Leica DMI5000 M

Inverted microscope for materials testing –
Operating Manual

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
Leica DMI5000 M

Inverted microscope
for materials testing

Operating Manual

Leica
MICROSYSTEMS

Copyrights

Copyrights

All rights to this documentation are held by Leica Microsystems Wetzlar GmbH. Reproduction of text or illustrations (in whole or in part) by print, photocopy, microfilm or other method (including electronic systems) is not allowed without express written permission from Leica Microsystems Wetzlar GmbH.

The term “Windows” may appear in the following text without further identification. It is, however, a registered trademark of Microsoft Corporation. The names of companies and products used herein may be trademarks of their respective owners.

The instructions contained in the following documentation reflect state-of-the-art technology and knowledge standards. We have compiled the texts and illustrations as accurately as possible. Nevertheless, no liability of any kind may be assumed for the accuracy of this manual's contents. Still, we are always grateful for comments and suggestions regarding potential mistakes within this documentation.

The information in this manual is subject to modification at any time and without notification.

Contents

1. Important Notes about this Manual	7	6.12 Installation and replacement of transmitted-light lamps: lamp housing 107 or 107/2	42
2. Intended Purpose of the Microscope ...	8	6.13 Installation of Condensers	44
3. Safety Notes	9	6.14 Installation of Eyepieces	47
3.1 General Safety Notes	9	6.15 Optional Accessories	48
3.2 Electrical Safety	9	6.16 Connection to the CTR6000/6500 electronics box	49
4. Overview of the Leica DMI5000 M	11	6.17 Connection to the CTR4000 electronics box	50
4.1 Specifications	11	6.18 Connection to the Computer	50
4.2 Glossary	14	6.19 Connection to the Power Supply	50
5. Unpacking the Microscope	18	7. Startup	51
6. Assembling the Microscope	20	7.1 Functional Principle	51
6.1 Assembly Tools	20	7.2 Switching on the Microscope	55
6.2 Installation of the Lamp Mount, Mirror Housing, Lamp Housing, Illumination Telescope	21	7.3 The LeicaScreen	56
6.3 Installation and Replacement of Incident-Light Lamps	23	7.4 The Function Buttons on the Stand	57
6.3.1 12 V 100 W Halogen Lamp	23	7.5 The SmartMove Remote Control Module	60
6.3.2 Hg and Xe Lamps	25	7.6 Incident Light	61
6.4 Installation of the DIC Module and DIC Objective Prisms	28	7.7 Adjusting the Light Sources	62
6.5 Installation of Objectives	30	8. Operation	65
6.6 Equipping the Reflector Turret	31	8.1 Switching on	65
6.7 Installation of the Incident-Light Polarizers and Analyzers	33	8.2. Contrast Methods	66
6.8 Installation of Stages	36	8.2.1 Bright Field (RL/TL)	66
6.9 Installation the transmitted-light illumination carrier	40	8.2.2 Dark Field (RL/TL)	67
6.10 Installation of Filters in the Illumination Arm	41	8.2.3 Polarization (RL/TL)	68
6.11 Installation of the Transmitted-Light Lamp Housing	41	8.2.4 Differential Interference Contrast (RL/TL)	69
		8.2.5 Fluorescence	70
		8.2.6 Phase Contrast (TL)	70
		8.3 Focusing	71
		8.4 Tubes	73

Contents

8.5	Optical Equipment	74
8.6	Graticules for Eyepieces	78
8.7	Eyepieces	81
8.8	Objectives	81
8.9	Stages and Specimen Displacement	84
8.10	Magnification Changer	85
8.11	Light Sources	85
8.12	Aperture and Field diaphragm	86
9.	Troubleshooting	87
10.	Care of the Microscope	91
10.1	Dust Cover	91
10.2	Cleaning	91
10.3	Handling Acids and Bases	92
11.	Essential Wear and Spare Parts	93
12.	Dimensions	94
13.	Abbreviations and Pictograms	95
14.	Index	97
15.	EU Declaration of Conformity	100

1. Important Notes about this Manual

Caution!

This operating manual is an essential component of the microscope, and must be read carefully before the microscope is assembled, put into operation or used.

This operating manual contains important instructions and information for the operational safety and maintenance of the microscope and accessories. It must therefore be kept safely for future reference.

A separate manual is available on CD-ROM covering the operation of the Leica Application Suite (LAS).

Text symbols, pictograms and their meanings:

(1.2)

Numbers in parentheses, such as "(1.2)", correspond to illustrations (in the example, Fig. 1, Item 2).

→ p. 20

Numbers marked with an arrow, such as → p. 20, refer to specific pages in this manual.

Caution!

Special safety instructions within this manual are indicated with the triangle symbol shown here, and have a gray background.

Caution! The microscope and accessories can be damaged when operated incorrectly.

Explanatory note.

*

Item not contained in all configurations.

2. Intended Purpose of the Microscope

2. Intended Purpose of the Microscope

The Leica DMI5000 M is designed for routine industrial and research applications.

The Leica DMI5000 M is the logical further development of Leica's proven inverted research microscopes. It is designed for use in metallographic laboratories for materials testing and research involving opaque and transparent structural materials.

The Leica DMI5000 M is suitable for universal deployment. All contrast methods such as dark field, bright field, phase contrast, DIC, polarization and fluorescence are integral to the microscope and can be adapted or changed quickly and easily. Variable illumination and imaging beam paths, as well as HCS optics, modular accessories and a comprehensive range of peripherals complement the Leica DMI5000M inverted research stand.

The above-named microscope complies with the Council Directive 98/79/EEC concerning in vitro diagnostics. They also conform to the Council Directives 73/23/EEC concerning electrical apparatus and 89/336/EEC concerning electromagnetic compatibility for use in an industrial environment.

Caution!

The manufacturer assumes no liability for damage caused by, or any risks arising from using the microscopes for other purposes than those for which they are intended or not using them within the specifications of Leica Microsystems Wetzlar GmbH. In such cases the declaration of conformity shall cease to be valid.

Caution!

This (IVD) device is not intended for use in the patient environment defined by DIN VDE 0100-710. Neither is it intended for combining with medical instruments according to EN 60601-1. If a microscope is electrically connected to a medical instrument according to EN 60601-1, the requirements defined in EN 60601-1-1 shall apply.

3. Safety Notes

3.1 General Safety Notes

This safety class 1 device is constructed and tested in accordance with
 EN 61010-2-101:2002,
 EN 61010-1:2001,
 IEC 1010-1:2001,
 Safety regulations for electrical measuring, control, and laboratory devices.

Caution!

In order to maintain this condition and to ensure safe operation, the user must follow the instructions and warnings contained in this operating manual.

Caution!

The devices and accessories described in this operating manual have been tested for safety and potential hazards.
 The responsible Leica affiliate or the main plant in Wetzlar must be consulted whenever the device is altered, modified or used in conjunction with non-Leica components that are outside of the scope of this manual.

Unauthorized alterations to the device or noncompliant use shall void all rights to any warranty claims!

3.2 Electrical Safety

General specifications

Leica CTR4000, CTR6000, CTR6500 electronics boxes

For indoor use only.
 Supply voltage: 90–250 V~
 Frequency: 50–60 Hz
 Power input: max. 290VA
 Fuses: T6.3 A
 (IEC 60127-2/3)
 Ambient temperature: 15–35°C
 Relative humidity: max. 80% to 30°C
 Overvoltage category: II
 Pollution degree: 2

Microscope

For indoor use only.
 Supply voltage: 90–250 V~
 Frequency: 50–60 Hz
 Power input: → CTRxxx
 Fuses: → CTRxxx
 Ambient temperature: 15–35°C
 Relative humidity: max. 80% to 30°C
 Overvoltage category: II
 Pollution degree: 2

ebq 100 supply unit*

For indoor use only.
 Supply voltage: 90–250 V~
 Frequency: 50–60 Hz
 Power input: max. 155VA
 Fuses: 2xT2A (IEC 127)
 Ambient temperature: 10–36°C
 Relative humidity: max. 80% to 30°C
 Overvoltage category: II
 Pollution degree: 2
 (see enclosed manual)

3. Safety Notes

Caution!

Power plugs may only be plugged into an outlet equipped with a grounding contact.

Do not interfere with the grounding function by using an extension cord without a ground wire. Any interruption of the ground wire inside or outside of the device, or release of the ground wire connection, can cause the device to become hazardous. Intentional ground interruption is not permitted!

Caution!

Peripheral devices with their own or separate power supplies that are connected to the microscope can have the same protective conductor potential by connecting them to the ground screw on the back of the electronics box. For connections without a ground connector, Leica Service must be consulted.

Caution!

Never use any fuses as replacements other than those of the types and the current ratings listed here. Using patched fuses or bridging the fuse holder is not permitted. The use of incorrect fuses may result in a fire hazard.

Caution!

The microscope's electrical accessory components are not protected against water. Water can cause electric shock.

Caution!

Protect the microscope from excessive temperature fluctuations. Such fluctuations can lead to the accumulation of condensation, which can damage the electrical and optical components.
Ambient temperature: 15–35°C.

Caution!

Before exchanging the fuses or lamps, be absolutely certain to switch off the main power switch and remove the power cable.

4. Overview of the Leica DMI5000 M

4.1 Specifications

Contrast Methods	<ul style="list-style-type: none"> incident light (IL): BF, DF, POL, DIC (ICR), Fluo transmitted light (TL): BF, DF, PH, POL, ICT
Incident light axis	<ul style="list-style-type: none"> automatic Illumination Manager (aperture, field diaphragm, intensity, IL/TL switching) automatic constant-color intensity control motorized shutter (switching speed < 50ms) lamp housing mount for up to 3 interchangeable light sources motorized 4-place filter turret
Transmitted Light Axis	<ul style="list-style-type: none"> automatic Illumination Manager (aperture, intensity, IL/TL switching) manual field diaphragm, manual filter selection manual shutter lamp housing mount for interchangeable lamp housings automatic, electronic condenser identification
Tube	<ul style="list-style-type: none"> ergonomic with or without camera port at left 2 switching positions: 100% VIS or 100% CAM eye spacing adjustment height and angle adjustment (30°–45°)
Magnification Changer	<ul style="list-style-type: none"> motorized 3 switching positions (choice of magnifications: 1x; 1.5x; 1.6x or 2.0x) effective on all camera ports
Objective Turret	<ul style="list-style-type: none"> motorized and coded 5x for objectives with M32 thread and 45 mm parfocal distance for DIC (ICR): motorized/coded Wollaston prism turret anti-vibration locking

4. Overview of the Instrument

Stages	<p>Fixed stage</p> <ul style="list-style-type: none">• ceramic-coated stage plate (248 mm x 204 mm) <p>manual 3-plate cross-stage</p> <ul style="list-style-type: none">• positioning range: 83 mm x 127 mm• optional inserts for a variety of applications, size of inserts: 160 mm x 110 mm (compatible with scanning stages) <p>Scanning stage IM 120 x 100 (motors on bottom)</p> <ul style="list-style-type: none">• 1 mm, 2 mm, 4 mm spindle pitch (higher resolution vs. higher speed)• optional inserts for a variety of applications, size of inserts: 160 mm x 110 mm
Condensers	<ul style="list-style-type: none">• motorized and coded• motorized aperture diaphragm• contrast methods: BF, DF, PH, ICT, Pol• automatic contrasting mode switching• condenser turret with 7 positions for contrast methods• condenser housing for S1-S28• condenser heads: S1/1.4 oil, S1/0.9 dry, S23/0.53, S28/0.55• condenser heads can be swung out• all condensers suitable for magnifications from 1.25x to 100x• with or without motorized or manual polarizer• with or without motorized or coded Wollaston prism disk
Z focus	<ul style="list-style-type: none">• mechanical or motorized (both coded)• 9 mm travel (1 mm above, 8 mm below the stage)• maximum travel speed: 5 mm/s• 5 focus steps: 0.05 μm; 0.1 μm; 0.7 μm; 1.5 μm; 5.0 μm for motorized version• electronic focus repositioning• automatic lowering prior to objective change• electronic parfocality
Observation ports	<ul style="list-style-type: none">• motorized and coded<ul style="list-style-type: none">• left side ports (100%, 80% or 50% transmission)• right side ports (100%, 80% or 50% transmission) <p>optional</p> <ul style="list-style-type: none">• top port with 2 switching positions<ul style="list-style-type: none">• 100% to eyepieces/ 0% to camera port• 0% to eyepieces/ 100% to camera port

4. Overview of the Instrument

Controls	<ul style="list-style-type: none">• 7 fixed control buttons for illumination and apertures• 7 variable function buttons behind the focus controls• 3 fixed control buttons for focus stops (motorized focus)• 2 focus hand wheels• 4 buttons for reflector cubes and shutters• 4 buttons for magnification changer and ports• SmartMove: ergonomic remote control module for x,y,z control and four additional variable function buttons																														
Electronics box Leica CTR4000 Leica CTR6000 Leica CTR6500	<ul style="list-style-type: none">• separate control unit for all motorized and electronic elements of the microscope such as:<ul style="list-style-type: none">• objective turret• focus• ports• magnification changer• reflectors• condenser• motorized stageswith<ul style="list-style-type: none">• power supply for 100 W halogen lamps• power supply for SmartMove <table><tr><td></td><td colspan="2">Z-drive</td><td colspan="3">stage</td></tr><tr><td></td><td>manual</td><td>motorized</td><td>manual</td><td>rack</td><td>scanning</td></tr><tr><td>CTR4000 Order No. 11 888 164</td><td>x</td><td></td><td>x</td><td></td><td></td></tr><tr><td>CTR6000 Order No. 11 888 109</td><td></td><td>x</td><td></td><td>x</td><td></td></tr><tr><td>CTR6500 Order No. 11 505 220</td><td>x</td><td>x</td><td></td><td></td><td>x</td></tr></table>		Z-drive		stage				manual	motorized	manual	rack	scanning	CTR4000 Order No. 11 888 164	x		x			CTR6000 Order No. 11 888 109		x		x		CTR6500 Order No. 11 505 220	x	x			x
	Z-drive		stage																												
	manual	motorized	manual	rack	scanning																										
CTR4000 Order No. 11 888 164	x		x																												
CTR6000 Order No. 11 888 109		x		x																											
CTR6500 Order No. 11 505 220	x	x			x																										
Interfaces	<ul style="list-style-type: none">• 2 x RS232C• 2 x USB• 4 x external/internal peripherals																														
Software tools	<ul style="list-style-type: none">• Leica Application Suite (LAS) for Windows™ 2000, XP with plug-ins for:<ul style="list-style-type: none">• microscope and camera configuration• microscope and camera control• image acquisition																														

4. Overview of the Instrument

4.2 Glossary

The Stand

Two basic Leica DMI5000 M stand types are available that can be combined into a wide range of microscope versions.

The basic modules of the Leica DMI5000 M stand are:

- DMI5000 M electronic stand
- integrated incident-light axis with motorized filter cube changer (4 positions)
- lateral photo ports, 100%, 80% or 50%
- optional top port at the left side of tube
- integrated motorized magnification changer

The individual versions and their components, differences and application areas are covered in this manual. All microscopy processes and the required Leica DMI5000 M accessories are explained in detail in the operation section of this manual.

Tube

The tube and tube lens are responsible for the primary image in conjunction with the objective. The tubes are integral parts of the stand, consisting of a basic body and a binocular section. The trinocular tube features an additional photo/TV port. A switchable mirror diverts 100 % of the light to the eyepieces or the photo port.

Eyepieces

The eyepieces provide an enlarged, virtual image of the intermediate image created by the objective. The eyepiece acts as a magnifying glass in this respect.

Intensity controls

The Leica CTR box contains a 12V 100W transformer for continuous brightness adjustment via the intensity controller.

The intensity can be adjusted using the controls described in 1.5, p. 16.

Focus wheel

The focus wheel permits the microscope image to be focused quickly and precisely. Focusing is realized via a vertical motion of the objective turret, with a stroke of 9 mm.

Incident-light reflector turret

The stand features an integrated incident-light axis with a motorized 4-position reflector holder.

Aperture diaphragm

The aperture diaphragm affects the resolution, depth of field and contrast of the microscope image.

Caution:

The aperture diaphragm in the illumination beam path is not intended for adjusting the brightness of the image. Use only the intensity controls or neutral optical attenuation filters for this purpose.

Condenser

The condenser is a lens system that collects light and directs it onto the specimen from above. The condenser ensures the utilization of the numerical aperture of the objective.

4. Overview of the Instrument

Condenser height adjustment

The markings on the transmitted-light illumination column indicate the heights to be set for specific condensers.

Stages and accessories

The stage is designed to hold the specimens to be observed. Manual stages, motorized 3-plate-stages and scanning stages are available for the Leica DMI5000 M.

Motorized objective nosepiece and objectives

The motorized objective nosepiece accommodates the microscope's objectives. The L objectives with their long working distances are designed especially for the varying thicknesses of the specimen inserts.

Microscope objective from 1.25 :1 to 100:1 magnification may be used. All Leica M32 objectives are compatible; M25 may be used with an adapter ring. For the specifications of Leica objectives, please refer to the current objective lists which are available from your Leica representative.

Incident-light axis

The incident-light axis unit consists of the illumination carrier (lamp, aperture diaphragm, field diaphragm) and the four-position reflector turret. A constant-color intensity control unit (CCIC) is always integrated.

Transmitted-light axis

The transmitted-light axis consists of an illumination arm and the transmitted-light illumination column. While the transmitted-light illumination column does not have a lamp housing (lamp housing is ordered separately and is mounted on

the column), it features a filter module for two manual swing-in filters and a manual switchable field diaphragm.

Filter

The filters are generally used to enhance the contrast of the specimen. They are installed in the illumination column. A variety of filters can be installed with ease.

Field diaphragm

The field diaphragm is used to realize Koehler illumination.

Lamp housing for transmitted light

The Leica DMI5000 M is available with the 107/2 and 107 lamp housings (both for 12 V 100 W halogen). For their descriptions and areas of deployment, please refer to the operation section of this manual. The letter L signifies a lamp housing designed for left-hand use.

Lamp housing for incident light

The Leica DMI5000 M is available with the 107 L (halogen) lamp housing. For their descriptions and applications, please refer to the operation section of this manual. The letter L signifies a lamp housing designed for left-hand use.

Leica CTR4000, CTR6000, CTR6500 electronics box

The Leica electronics boxes contain the power supply for the lamp and electronics boards required to control the motorized functions of the stand.

4. Overview of the Instrument



Fig. 1 Leica DMI5000 M left view

- 1 Eyepiece
- 2 Eyepiece tube
- 3 Top port
- 4 LeicaDisplay
- 5 Light intensity
- 6 Field diaphragm
- 7 TL/IL switching
- 8 Aperture diaphragm
- 9 Focus wheel
- 10 Variable function buttons
- 11 Right side port
- 12 Lamp mount with incident-light lamp housing
- 13 Filter slide, e.g. for ND4
- 14 Polarizer mount
- 15 DIC objective prism disk

Fig. 2 Leica DMI5000 M front view

- 1 LeicaDisplay (two-color, high contrast)
- 2 Front control panel
- 3 Top-port switching (manual)
- 4 Top port
- 5 Motorized port switching (optional)
- 6 Shutter control
- 7 Magnification changer switching
- 8 Reflector switching

4. Overview of the Instrument


Fig. 3 Leica DMI5000 M right view

- 1 Incident light lamp housing
- 2 Aperture diaphragm centering
- 3 Field diaphragm centering
- 4 Centering window (lamp adjustment)
- 5 Analyzer slot
- 6 Right side port
- 7 Variable function buttons
- 8 Focus wheel
- 9 Opener for reflector drawer
- 10 Reflector drawer
- 11 Objective turret with objectives
- 12 Stage with attachable mechanical inserts
- 13 Focus buttons (optional - motorized version)
- 14 Transmitted illumination (optional)


Fig. 4 SmartMove remote control module

- 1 Move stage in x
- 2 Move stage in y
- 3 Focus
- 4 Variable function buttons (preassigned at factory)

5. Unpacking the Microscope

The microscope is delivered in several packages.

The **stand package** contains the following components:

- Stand with integrated incident-light axis, objective turret and tube
- Illumination arm
- Specimen stage
- CD with Leica Application Suite (LAS) software package
- Instructions and list of microscope presets (identification sheet)

The **system package** contains the microscope's accessories:

- Eyepieces
- Objectives
- Condenser
- Lamp housings with accessories
- Assembly tools
- Additional accessories such as reflector cubes, etc. depending on feature set

The electronics box, the SmartMove remote control module and the external ebq 100 supply unit, where relevant, are provided in separate packages.

Please carefully compare the contents of the delivery to the packing slip, delivery note or invoice. We urgently recommend storing a copy of these documents with the manual to ensure that you have information on the time and scope of delivery handy for subsequent orders or service work. Please ensure that no small parts remain in the packing material. Parts of our packing material are marked by symbols to simplify recycling.

First, carefully remove all components from the transportation and packaging materials.

Caution!

Do not put the instrument into operation in the event of visible damage to the components or packing material.

Note:

If at all possible, avoid touching the lens surfaces of the objectives. If fingerprints do appear on the glass surfaces, remove them with a soft leather or linen cloth. Even small traces of finger perspiration can damage the surfaces in a short time. See the chapter "Care of the Microscope", for additional instructions → p. 91.

5. Unpacking the Microscope

Caution!

Do not connect the microscope or peripherals to an AC power source at this time under any circumstances!

Installation location

Work with the microscope should be performed in a dust-free room, which is free of oil vapors and other chemical vapors, as well as extreme humidity. At the workplace, large temperature fluctuations, direct sunlight and vibrations should be avoided. These may adversely affect measurements and long-term observations.

Allowable ambient conditions

Temperature 15–35°C

Relative humidity maximum 80% up to 30°C

Microscopes in warm and warm-damp climatic zones require special care in order to prevent the build up of fungus.

See the chapter “Care of the Microscope” → p. 91, for additional instructions.

Caution!

Electrical components must be placed at least 10 cm from the wall and away from flammable substances.

Transport

For shipping or transporting the microscope and its accessory components, the original packaging should be used.

As a precaution to prevent damage from vibrations, the following components should be disassembled and packaged separately:

- Unscrew the objectives.
- Remove the eyepieces.
- Remove the condenser if present.
- Remove the specimen stage.
- Remove the transmitted-light arm if present.
- Remove the lamp housings.
- Remove the lamp housing mount.
- Disassemble the burner of 106 z lamp housing.
- Remove the reflector cubes.
- Remove all moving or loose parts.

6. Assembly

6. Assembling the Microscope

The microscope components* are logically assembled in this order:

- Lamp housing mount, mirror housing, lamp housing
- Transmitted-light illumination carrier*
- Incident-light lamps
- Transmitted-light lamps*
- DIC module and DIC objective prisms*
- Objectives
- Fixation of reflector cubes onto the turret
- Polarizer and analyzer*
- Specimen stage
- Condenser with condenser head*
- Eyepieces

6.1 Assembly Tools

If possible, the microscope should be assembled and set up with the assistance of Leica sales or service personnel.

A small number of universal screwdrivers which are included in the scope of delivery are required for assembly (Fig. 5).

Fig. 5 Assembly tools

- 1 Phillips screwdriver*
- 2 3 mm Allen key
- 3 1.5 mm centering key*
- 4 2 mm centering key*
- 5 3 mm hex key*
- 6 2.5 mm hex key* (short type)
- 7 2.5 mm hex key*

* depending on scope of delivery

6. Assembly

6.2 Installation of the Lamp Mount, Mirror Housing, Lamp Housing, Illumination Telescope

Caution!

Light sources pose a potential irradiation risk (glare, UV-radiation, IR-radiation). Therefore, lamps have to be operated in closed housings.

Ensure that the lamp housing has been disconnected from the power supply. Unplug the power plug and the power supply during assembly.

During assembly work on xenon burners, always wear the supplied protective gloves and face protection (Fig. 6) (risk of explosion).

Never touch the glass parts of the burner with bare hands.
Never look directly into the beam path (risk of damage to the eyes).

Fig. 6
Protective gloves and mask

- Insert the lamp mount or mirror housing in the rear panel and fasten with socket screws.

Fig. 7 Leica DMI5000 M rear view

- 1 Installation point for lamp housing mount or mirror housing
- 2 Threaded holes for mounting screws of the lamp housing mount or mirror housing




Fig. 8 Lamp mount
1 Lateral lamp housing mount
2 Mounting screws (socket screws)

6. Assembly

- Place lamp housings 107/2, 107 L and 106 z on the lamp housing mount and tighten mounting screw. (Fig. 9).
- We recommend the use of an illumination telescope when using gas discharge lamps in conjunction with fluorescence methods. It is inserted between the lamp mount and the lamp housing 106 z and enlarges the image of the lamp's focal point by a factor of 2x in the entrance pupil of the objective. This results in significantly higher illumination intensity when working with fluorescence.

Fig. 9 Rear view Leica DMI5000 M with lamp housing mount

6.3 Installation and replacement of incident-light lamps

Caution!

Make sure to follow the safety notes on page 21.

6.3.1 Changing the 12 V 100 W halogen lamp

- Disconnect the lamp and lamp housing from the power supply.
- Unplug the power plug and power supply.

Lamp housing 107 L

- Remove the fastener screw (10.4) on the housing and lift the housing off.
- Move the collector (10.2) forward and pull the defective 12 V 100 W lamp out of its socket. (11.1)

Caution: surface may be hot!

- Press the new lamp and its dust cover home into the socket.

Attention: Do not remove the new lamp's dust cover until you have installed the lamp. Avoid fingerprints on the lamp and remove them if necessary.

- Close the lamp housing.

Fig. 10 Lamp housing 107 L

- 1 Mounting dovetail ring
- 2 Collector adjustment
- 3 Horizontal and vertical lamp adjustment
- 4 Locking screws for housing



Fig. 11 Lamp housing 107 L opened

- 1 Mount with halogen bulb
- 2 Collector
- 3 Threaded hole for housing

6. Assembly

Lamp housing 106 z L

- Release the mounting screws of the cover. (12.4, 9)
- Loosen the contact plug in its socket (12.11) somewhat and flip the cover upwards. (12.1)
- Move the collector forward and pull the defective lamp upwards out of its socket. (12.2, 12.3, 14.1)
- For added convenience, the lamp mount can be removed from the lamp housing. Remove the mounting screws from the lamp mount (12.10) and pull the lamp mount out.

Caution: surface may be hot!

- Press the new lamp and its dust cover home into the pin base.
Attention: Do not remove the new lamp's dust cover until you have installed the lamp. Avoid fingerprints on the lamp and remove them if necessary.
- Close the lamp housing.

Fig. 12 Lamp housing 106 z L

- 1** Cover raised
- 2** Collector
- 3** 12 V 100 W lamp or gas discharge lamp in mount
- 4, 9** Cover fixation
- 5** Reflector
- 6, 8** Adjusting screws for x-y reflector
- 10** Locking screws/lamp mount
- 11** Socket for contact plug


Fig. 13 12 V 100 W lamp mount with halogen bulb

6. Assembly

6.3.2 Installation and replacement of Hg and Xe lamps

Lamp housing 106 z

This lamp housing is suitable for use with a 12 V 100 W halogen lamp or a variety of gas discharge lamps.

Caution!

Make sure to follow the instructions and safety notes of the lamp supplier.
Before changing lamps allow at least 30 mins for cooling down!

Inserting gas discharge lamps (Hg and Xe) in the 106 z lamp housing

Hg and Xe lamps are powered by separate supply units.

Please also read the separate instruction manual provided with these supply units.


The following gas discharge lamps may be used and require different supply units and lamp mounts (Fig. 15, 16):

Type	Typical bulb life*
100 W high-pressure mercury burner (direct current)	200 hrs.
100 W high-pressure mercury burner (direct current, type 103 W/2)	300 hrs.
75 W high-pressure xenon burner (direct current)	400 hrs.

* Please observe the data sheets of the lamp manufacturer.

Fig. 14 Lamp housing 106 z L with Hg 100 W lamp

- 1 Collector focusing
- 2 Vertical lamp adjustment
- 3 Horizontal lamp adjustment
- 4 Hg lamp mount
- 5 Reflector adjustment (not visible)

6. Assembly

Caution!

Make sure to follow the safety notes on page 21.

- To open the 106 z lamp housing, unscrew the fastening screws on the cover (18.8). Loosen the contact plug somewhat and pull it out of the socket (18.9). Flip the cover up (18.1).
- Loosen the mounting screws (18.8) on the lamp socket and pull the socket out.
- Remove the transport anchorage (red plastic rod in place of the burner) in the lamp mount. To do so, remove the lower clamp (15.1, 16.1). Pull up the cooling element (16.3, 15.3) and turn it to the side. Detach the lower clamp system (16.2, 15.2) and remove the transport anchorage.

Caution!

Do not remove the burner's dust cover until you have installed the lamp. Avoid fingerprints on the lamp. Sweat from your fingers on the glass will shorten the life of the lamp significantly.

- Install the burner in reverse order.

Caution!

Xe 75 burner:

Remove the burner's dust cover (16.4) after you have installed the burner.

Fig. 15 Lamp mounts for Hg 100 gas discharge lamp

- 1 Upper clamping system
- 2 Lower clamping system
- 3 Cooling element


Fig. 16 Lamp mounts for gas discharge lamp Xe 75

- 1 Upper clamping system
- 2 Lower clamping system
- 3 Cooling element
- 4 Protective cover of Xe 75 burner

6. Assembly

- Insert the lamp mount, with the burner installed, into the lamp housing and tighten it with the screws (18.8).
- Connect the lamp housing to the external power supply (17.1).
- Test the adjustment of the collector (18.2):
Do not touch the power supply while performing these actions. When closing the lamp housing, ensure that the pins of the contact plug engage in their sockets (18.9). Tighten the screws of the cover and press the contact plug home.
- Place the lamp housing in the incident light lamp housing mount (Fig. 9) and fasten it with the clamping screw on the side.

Caution!



The burner must be adjusted immediately after lighting.

Fig. 18 106 z lamp housing (on the side, open)

- 1 Cover raised
- 2 Collector
- 3 12 V 100 W lamp or gas discharge lamp in mount
- 4 Reflector (mirror)
- 5, 6, 7 Adjusting screw for x-y reflector
- 8 Locking screws for lamp mount
- 9 Socket for contact plug

Fig. 17 Rear panel of ebq 100 supply unit

- 1 Lamp connection

6. Assembly

6.4 Installation of the DIC Module and DIC Objective Prisms

If your microscope is not equipped with DIC, please continue with Chapter 6.5, p. 30.

In the Leica DMI5000 M microscope, the DIC prisms are already installed in the DIC disk below the objective turret (Fig. 20).

Proceed as follows when making changes to the IC prism disk:

- Remove the front cover (Fig. 19, → 1.15, p. 16) under the objective turret after releasing the socket screws.

- Insert the DIC prism disk in its receptacle (Fig. 21) and tighten the two socket screws (Fig. 22).

Note: insert the prism disk with the electronics board facing down.

Fig. 19 Front cover, DIC prism disk

Fig. 21 Inserting the DIC module


Fig. 22 Mounting the DIC module

Fig. 20 DIC objective prism turret (coded and motorized)

- 1 IC objective prism in frame
- 2 Identification letter (ID)
- 3 Orientation pin

6. Assembly

Replacing individual IC prisms:

- Release the two socket screws and remove the prism disk.
- Place the prism against the stop pin, place the washer between the screw and the prism, and tighten gently to prevent undue tension. Insert the prism so that its identifying letter, e.g. D, is facing upward and is legible.
- After installing the prisms, replace the prism disk in its receptacle.

Fig. 23 IC objective prism

- 1** Objective prism in frame
2 Screw and washer

6. Assembly

6.5 Installation of objectives

BD (bright and dark-field) objectives are screwed directly into the M32 x 0.75 mm threaded openings of the objective turret (Fig. 25), clockwise from the lowest to the highest magnification.

Use 32/25 or 32/RMS spacer rings for M25x0.75 mm or RMS-threaded objectives. (Fig. 24).

The positions in the objective turret disk are numbered (Fig. 25). Depending on your equipment, the individual objectives have already been assigned to specific positions at the factory.

For details on the exact positions of the objectives, please refer to the enclosed identification sheet.

! Caution:

Close vacant threads in the nosepiece with dust protection caps!

Please note that the front lenses of the objectives point upward and are therefore more vulnerable to contamination than those of upright microscopes.

Check the front lenses for cleanliness frequently.

Note:

We recommend running a parfocality compensation via the Leica Application Suite (LAS) software for stands with motor focus.

Fig. 24 Spacer rings


Fig. 25 Objective turret, not loaded

Fig. 26 Objective turret, loaded

6. Assembly

6.6 Equipping the Reflector Turret

The 4-fold reflector turret drawer is located on the right side of the stand. Before opening this drawer, remove the cap below the drawer covering the analyzer slot. Remove the analyzer if it is already in the slot.

The replacement of individual cubes is more convenient with the microscope switched on. Turn the position to be replaced outward and place the cube in the correct mount.

The reflectors may also be installed immediately, with the microscope switched off. Press the white button next to the drawer. The drawer will glide out into its initial position. The inner disc for the reflector cubes can be turned in this position to access the desired position.

Up to 4 reflectors can be installed in the disk.

Fig. 29 Opening the 4-fold reflector turret drawer

Fig. 30 Open reflector drawer

Fig. 31 Inserting or removing a reflector cube

Fig. 27 Reflector cube, front side

Fig. 28 Reflector cube, back side

6. Assembly

The positions in the reflector turret are numbered. Positions 1 and 4 are equipped with quick releases. Reflectors can be snapped into these positions as described below. Positions 2 and 3 use screw fittings.

Depending on your equipment, the individual reflector cubes have already been assigned to specific positions at the factory. For details, check the identification sheet included with your order.

Now open the drawer until it clicks into its end position. The disk can be locked into this position to prevent it from turning further.

You may now insert a reflector cube. Proceed as follows:

- With the holder facing you squarely, insert the reflector cube into the holder in accordance with the included identification sheet.

- The reflector cubes are suitable for both upright and inverted microscopes. When using them with inverted microscopes, insert them so that the writing is upside down along the lower edge.

To do so, place the reflector cube on the **left** side and press it to the **right** into the mounting (Fig. 31).

- Ensure that the cube is correctly seated. A loose cube can block the disk or be destroyed by the turning disk.
- For the next cube, unlock the disk so that it is once again free to turn. Once you have reached the next position, relock the disk. Continue in this way for all of the cubes.
- Once all filter and reflector cubes have been inserted, close the drawer and replace the analyzer or cap.

Fig. 32

Reflector BF
1 Neutral plane
face glass
divider
2 Neutral
gray filter N16

Fig. 33

ICR reflector

Fig. 34

DF reflector
(threaded)


Fig. 35

Smith reflector
(threaded)

6.7 Installation of the Incident-Light Polarizers and Analyzers

Incident-light polarizer R/P in slider (Fig. 36)

The vibration axis of the polarizer can be changed as required,

- 0° = East-West
- 45° = Diagonal
- 90° = North-South

The vibration axis 0° = East-West is intended for use with the DMI5000 M.

Pull the polarizer holder (36.1) out of the slider and insert it in the appropriate detent.

In conjunction with analyzer IC/P at 90° = North-South, the polarizer is now in cross position.

Further incident-light polarizers in sliders:

Incident-light polarizer, rotatable 90° with adjustable lambda plate (Fig. 37)

The vibration axis 0° = East-West is intended for use with the DMI5000 M.

Preselect the vibration axis 0° = East-West at the selector disk. (37.2)

In conjunction with analyzer IC/P at 90° = North-South, the polarizer and analyzer are now in cross position.

Polarization color contrast ensues when polarizers are crossed and the lambda plate is turned using the knurled button (37.3).

The rotating range of the lambda plate is 14°. Especially advantageous for color contrasting of anisotropic metal surfaces.

Fig. 36 Incident-light polarizer R/P in slider

- 1 Polarizer holder, removable
- 2 Polarizer slider
- 3 0°, 45°, 90° detents on rear



Fig. 37 Incident-light polarizer, rotatable 90° with rotatable lambda plate

- 1 Polarizer slider
- 2 Polarizer adjuster disk
- 3 Knurled knob for adjusting lambda plate

6. Assembly

Incident-light polarizer L ICR/P, with lambda plate, reversible (Fig. 38)

For use with the DMI5000 M with fixed 0° = East-West vibration axis.

The lambda plate has a fixed 45° diagonal orientation and can be enabled or disabled by reversing the polarization slider by 180° . The polarization slider is faceted on both sides for reversal. (38.2)

Activating the lambda plate results in an interference color contrast with a fixed path difference of one light wavelength. (red first order)

The polarizer L ICR/P is also used for grayscale and color scales in incident-light interference contrast.

Insert the polarizer slider in the polarizer slot (1.14, p.16) so that the engraving points to the user.

The polarizer slider features two detent notches that ensure the correct position of the slider in the beam path, On or Off.

Caution:

When using highly intense gas discharge lamps such as the Xe 75 W, you must insert the Pol protection filter 504079 in the lamp housing (39.1) to prevent the destruction of the polarizer foil.

Fig. 38 Incident-light polarizer L ICR/P with lambda plate

- 1 Polarizer slider
- 2 Facetted edges, both sides



Fig. 39 Inserting the Pol protection filter

- 1 Pol protection filter

6. Assembly

Incident-light analyzer

- Remove the cap (Fig. 40) on the right side of the stand (under the reflector drawer).
- Insert the analyzer into the receptacle until it latches in place (Fig. 41.1).

Analyzer IC/P (Fig. 43)

for polarization contrast and interference contrast. Vibration axis 90° North-South, rotatable from 0 to 180°; drum head graduation: 1 graduation mark = 5°.

The polarization is in cross position when combined with one of the polarizers mentioned above.

Insert the analyzer in its slot with the writing facing up (3.5, p. 17).

Fig. 40 Analyzer slot cap

Fig. 43 Analyzer IC/P, 180°


Fig. 41 Inserting the analyzer

- 1 Slot
- 2 Analyzer

Fig. 42 Inserting the analyzer

6. Assembly

6.8 Installation of Stages

A wide range of stages are available. The most important are the following:

- Fixed stage (248 mm x 204 mm):
- manual 3-plate cross-stage positioning range: 83 mm x 127 mm
- manual rotating stage
- Motorized 3-plate cross-stage positioning range: 83 mm x 127 mm
- scanning stage IM 120 x 100 (motors on bottom)

The assembly of these stages is identical. The stages are solidly attached to the microscope by three screws. In the case of a fixed stage, a mechanical stage can be attached at the right or left (Fig. 46). The attachable mechanical stage is supplied in a separate package.

Multiple-plate stages are supplied separately. Like the fixed stages, these stages are mounted as follows:

- If the screws for the stage are already in the stand, remove them first. In most cases, the screws will be found in the packing material of the stand.

Caution!

The screw lengths may vary. When using screws of different lengths, use the shorter of the three screws in the front hole and the equally long ones in the rear holes.

Fig. 44 Fixed stage (normal)


Fig. 45 Mechanical 3-plate stage


Fig. 46 Attachable mechanical stage for fixed stage

6. Assembly

- Use a clean cloth to remove dust and packing material residue from the stand's contact surface for the stage.
- Align the stage so that the pair of holes faces back toward the illumination axis and the single hole faces forward toward the tube.
- Align the mounting holes in the stage with the holes in the support surface. If the holes are covered in the case of 3-plate cross-stages or scanning stages, please shift the upper stage plate until the opening becomes visible.
- First, tighten the single front screw with the included 3 mm hex screwdriver. Be sure to use the shortest of the three screws in the front hole, as an excessively long screw can interfere with the focus travel. (If you have a rotating stage, please continue reading under "Rotating Stage and Insert Frame for Cover-slips").
- Next, firmly tighten the two rear screws.
- Finally, give the front screw a final firm tightening.

Fixed stage

Attachable mechanical stages designed to accept a variety of specimen holders are also available for fixed stages. They may be attached to the fixed stage on the left or right as required.

Two screws are located on the underside (left or right) of the attachable mechanical stage. Tighten these screws in the threaded holes on the underside of the fixed stage with the 3 mm hex screwdriver. Retighten these screws from time to time after frequent use.

The attachable mechanical stage has been preadjusted in the factory. In the event that the attachable mechanical stage runs out of focus when moving from right to left, this can be corrected by Leica's technical service.

Next, remove one of the ordered insert frames from their packaging and place the insert frame into the precise retention system. The stage, the attachable mechanical stage and the insert frame are now ready for use.

Some (not all) inserts are provided with self-adhesive scales to permit the coordinates to be read.

Apply these scales to the recesses of the attachable mechanical stage.

6. Assembly

Motorized 3-plate or scanning stages

Motorized 3-plate stage: after installing the stage, connect the included stage cable first to the socket on the stage, then to the Leica CRT6000/CTR6500 box. The correct place on the box is called "XY Stage". (→ Table p. 13)

A variety of inserts are available for the 3-plate and scanning stages. Install these inserts diagonally from above into the corner with the spring clips. The insert will click into place when seated properly.


Caution:

Press the spring clip into place only from the side.

Do not press the insert onto the spring clips diagonally from above, as the insert will not be aligned parallel to the stage and may be bent in the process.

Fig. 47 a, b Mounting screws for 3-plate cross-stage

6. Assembly

Rotating Stage and Insert Frames for Coverslips

The rotating stage (Fig. 48) is also mounted with 3 screws (48.2). Rotate the stage to make all of the threaded holes accessible. Insert the screws (48.2).

! Caution:

Use additional washers (48.3) for the rear holes. Tighten the screws only lightly, as the rotating stage must be centered first: insert the adjusting aid and the focusing telescope (Fig. 50) into the rotating stage for this purpose. Move the stage until the bright circle is in the middle of the field of view. Next, tighten the stage, replace the focusing telescope with an eyepiece and remove the adjusting aid.

To insert glass slides in insert frames (49.1), press on the center of the leaf spring (49.2) and insert the coverslip in the direction of the arrow. Clamp the insert frame in the attachable mechanical stage (48.1).

Fig. 48 Rotating stage

- 1 Attachable mechanical stage
- 2 Screws for stage mounting
- 3 Washers


Fig. 49

- 1 Insert frame for coverslips
- 2 Leaf springs


Fig. 50 Focusing telescope

6. Assembly

6.9 Installation of the Transmitted-light Illumination Arm

Wipe the installation surface (51.3) with a dry cloth. Tip the illumination arm (51.1) back slightly and install it so that the pin (51.2) engages the groove in the support surface (51.4).

Set the TL illumination arm upright and fasten it with the 4 screws.

When fastening the illumination arm, do not hold it so as to ensure its optimal alignment with the optical axis.

The tilt angle of the illumination arm can be varied with the knurled screw (52.1) or fixed vertically. Connect the electronics cable to one of the sockets, EXT1 - EXT4.

The transmitted-light lamp housing for 12 V 100 W halogen lamps is a separate component. For instructions on replacing the halogen lamp → Ch. 6.12, p. 42.

Fig. 52 Transmitted-light illumination arm, rear side

- 1 Knurled locking knob of the transmitted-light illumination arm
- 2 Connector cable for the electronics box





Fig. 51 Installing the transmitted-light illumination carrier

- 1 Transmitted-light illumination arm
- 2 Transmitted-light illumination arm pin
- 3 Support surface
- 4 Support surface groove
- 5 Support surface groove
- 6 EXT1-EXT4 sockets
- 7 Electronics cable

6.10 Installation of Filters in the Illumination Arm

The Leica DMI5000 M is generally equipped with a holder for two 40 mm dia. filters. The filters are installed at the factory. To change filters yourself, proceed as follows:

- Release the screw (53.1) and remove the cover.
- Place the filter in the holder.
- Place the cover on the transmitted-light illumination arm and fasten with the locking screw.
- Mark the levers with the provided adhesive labels.

Fig. 53 Unscrewing the filter holder cover and inserting filters in the transmitted-light illumination arm

1 Screw



Fig. 55 Lamp housing cabling (cable duct)

6.11 Installation the Transmitted-light Lamp Housing

- Place the lamp housing in the transmitted light lamp housing mount (Fig. 54) and fasten it with the clamping screw on the side.
- Thread the cable through the transmitted-light illumination arm (Fig. 55).
- Connect the lamp housing cable to the power supply for transmitted light on the Leica CTR4000, CTR6000 or CTR6500 electronics box (Fig. 56.1).

For instructions on changing the lamp, please see Chapter 6.12.

Fig. 54 Mounting the lamp housing on the transmitted-light illumination arm

Fig. 56 Connecting the lamp housing to the Leica electronics boxes

6. Assembly

These instructions also apply to installing an Hg lamp on the transmitted-light axis. For descriptions of the lamp housings and replacement of the burner, please see Chapter 6.3, → p. 25ff.

6.12 Installation and Replacement of the Transmitted-light Lamps: Lamp Housing 107 or 107/2

This lamp housing is used with a 12 V 100 W halogen lamp, which is already mounted. In case the lamp has to be removed:

Changing the 12 V 100 W halogen lamp

Caution!

Ensure that the lamp housing has been disconnected from the power supply. Unplug the power plug and the power supply during assembly.

Caution!

Light sources pose a potential irradiation risk (glare, UV-radiation, IR-radiation). Therefore, lamps have to be operated in closed housings.

- Remove the fastener screw on the housing (Fig. 57).

Fig. 57
Lamp housing 107/2
Releasing the
fastening screw

Fig. 58
Removing housing

6. Assembly

- Lift the housing off (Fig. 58).
- Remove the lamp.

Caution!

Do not remove the new lamp's dust cover until you have installed the lamp. Avoid fingerprints on the lamp.

- Insert the new 12 V 100 W lamp (Fig. 60) with the dust cover straight into the socket until it stops. Be sure that the lamp is inserted straight.
- Remove the lamp's dust cover.
- Replace the housing and fasten it in place using the fastening screw.

Fig. 59
Lamp housing 107/2
opened

- 1 Mount with halogen lamp
- 2 Collector




Fig. 60
Inserting
lamp with
cover
a right
b wrong

6. Assembly

6.13 Installation of Condensers

All condensers of the Leica DMI5000 microscope are equipped with a 7-position turret disk that can be equipped with light rings for dark-field (DF) or IC prisms for transmitted-light interference contrast (DIC).

Light rings and condenser prisms are generally already installed in the turret at the factory, making the following assembly steps unnecessary.

To retrofit IC prisms, follow the instructions below:

Installation of IC prisms

Ensure that the instrument is switched off before changing IC prisms. To retrofit IC prisms, start by removing the complete condenser from the transmitted-light axis. Simply release the socket screw on the right side of the condenser mount (Fig. 68) to quickly and easily remove the condenser thanks to Leica's condenser quick-release mount. Open the condenser cover (65.1) at the top right side. You now have access to the various numbered openings for the IC prisms.

- IC prisms can only be inserted in the condenser disk's large receptacles with guide grooves.
- Select a position and ensure that the two mounting screws have been released to the point that they no longer extend into the position. When inserting prisms, remove one screw fully, as a single screw is sufficient for centering. The counteracting pressures of both screws may destroy the prism.

Fig. 61 Condenser base S1-S28


Fig. 62
Condenser head S1

Fig. 63
Condenser head S28

Fig. 64 Condenser prisms

6. Assembly

- Next, take the special condenser tool (Fig. 66.1).
- If possible, install the prisms 0 to 3 in ascending order. The numbering can be found on the edge of the crown gear.
- Turn the right-hand centering screw back fully with the adjusting key (66.2). To prevent the condenser disk from turning further, insert the adjusting key (66.2) into the left-hand centering screw of the disk. It may protrude a maximum of 1 mm into the opening.
- Grasp the prism to be installed with the condenser tool (the lettering must face upward and be legible) so that the tab of the prism ring is positioned to the center of the tool's cam and the upper edge of the prism is lying flat in the holder of the tool. The numbers K2 to K16 should be positioned toward the end of the tool. Press the cheeks of the tool to grasp the prism (Fig. 67).
- Two guide hooks are located on the underside of the prisms. These must fit into the two grooves of the opening. Insert the prism (holding the condenser tool angled slightly upward and at a 90° angle to the housing) so that the mount fits under the spring clip of the retainer (Fig. 67).

! Caution:

Do not press the spring clip down under any circumstances. This can destroy the clip or result in an unstable position of the prism.

Fig. 66 Open condenser

- 1 Condenser tool
- 2 Adjusting key


Fig. 65 Condenser

- 1 Condenser cover
- 2 Centering opening


Fig. 67 Inserting the prism

The designation must be visible when installed and oriented toward the center of the condenser. DIC images are not possible otherwise.

6. Assembly

Turn the prism to ensure that it snaps into position and release the tool.
Remove fingerprints or dust from the prism with care.

- Use the left centering screw to roughly center the prism. The right centering screw must **not** restrict the range of adjustment under any circumstances.
- The ID, e.g. K10, must be visible when the prism is installed and must be oriented toward the center of the condenser. The prism can also be inserted with the ID facing outward, but DIC images will not be available with the prism in this position.
- Note the number of the opening and the prism designation for entry into the Leica Application Suite (LAS).
- Repeat the process for the remaining prisms.
- Remove the adjusting key and close the condenser.
- Perform the fine adjustment with the telescope after switching the unit on.


Installation of Condensers

The installation procedure is identical for all condensers, S1 to S28.

Release the socket head screw at the right side of the condenser holder. Place the condenser on the retaining pins of the illumination arm and move the condenser to the correct height. Use the markings on the column and condenser to determine the correct position.

Once you have reached the correct position, tighten the socket head screw.

Fig. 68 Installation of condenser on transmitted-light illumination arm

Condenser heads

Four different condenser heads are available:

- 1) S1/1.40 oil
- 2) S1/0.90 dry
- 3) S23/0.53
- 4) S28/0.55

Condenser heads 3 and 4 are screwed directly into the condenser body. A spacer ring (70.2) must be screwed into the thread at the bottom of the condenser body prior to installing condensers 1 and 2. The S1 condenser heads then fit into this ring.

The S70 condenser is delivered complete with a condenser head, making additional assembly unnecessary.

6.14 Installation of Eyepieces

The eyepieces are inserted into the eyepiece tubes.

Note:

We recommend running a teach-in via the Leica Application Suite (LAS) software when using eyepieces not included in the scope of delivery. This will ensure that the total magnification shown in the LeicaScreen is correct.

Fig. 69 Eyepieces

Fig. 70 Installation of condenser heads S1

- 1 Condenser base
- 2 Spacer ring
- 3 Condenser head


Fig. 71 Installation of condenser head S28

6. Assembly

6.15 Optional Accessories

Camera

Connecting a camera

A camera can be installed using a C-mount or Vario mount.

- Screw the camera onto the C-mount.
- Place the C-mount or Vario mount onto one of the camera ports and secure it with the locking screw at the side.

Note:

When using a C-mount or Vario mount, run a teach-in via the Leica Application Suite (LAS) software.

Connecting multiple cameras

Two or more cameras – for example a digital and an analog camera – can be adapted as required.

- When using a DC type camera, connect the camera to the PCI card of your PC.
- When using a DFC type camera, connect the camera to the FireWire card of your PC.

Note:

Please read the separate operating manual of your digital camera.

Fig. 72 C-mount 0.63x

Fig. 73 C-mount 0.5x

6.16 Connection to the CTR6000/6500 Electronics Box

Note:

The CTR6000/6500 electronics box must **not** be used with other stands. The serial number of the associated stand has been recorded on the back of the electronics box.

- Connect the **Microscope** (74.6) socket to the back of the stand (75.5) using the 25-pin microscope cable.

- Connect the SmartMove remote control module to the **XYZ-Control** socket (74.5).
- Connect the motorized stage, if present, to the **XY-Stage** socket (74.2).
- Connect the lamp power cable to the **12 V, max 100 W** socket (74.7).

Caution!

Ensure that the plugs are correctly inserted and secured to prevent overheating of the sockets.

Fig. 74 Rear view of Leica CTR6000

- 1 AC power socket
- 2 **XY Stage** socket for motorized stage
- 3 **Direct interface** socket optional
- 4 **Z Control** for separate focus control
- 5 **XYZ Control** for SmartMove
- 6 **Microscope** socket for microscope
- 7 **12 V, max 100 W** for the lamp power cable of stand
- 8 DL: reset button



Fig. 75 Rear view of stand

- 1 2 x RS232 ports
- 2 2 x USB
- 3 4 x EXT.
- 4 XYZ control for SmartMove
- 5 Connection to the Leica CTR6000/6500 electronics box

6. Assembly

6.17 Connection to the CTR4000 Electronics Box

- Connect the **Microscope** (76.2) socket to the back of the stand (75.5) using the 25-pin microscope cable.
- Connect the lamp power cable to the **12 V, max 100 W** socket (76.3).

6.18 Connection to the Computer

Note:

To start the Leica Application Suite (LAS), ensure that the COM1 serial port is not in use by another program or driver. This is frequently the case when using Palms or other PDAs or when using external modems or other devices. The

devices in question must therefore **always** be disabled before using the Leica Application Suite (LAS) software.

- Please use the included serial cable. Connect the COM1 port of your PC with the RS232C port (75.1) on the back of the stand.

6.19 Connection to the Power Supply

- Once all installation work is complete, connect the Leica electronics box to an AC power outlet with the included power cable (socket 74.1, 76.1).
- If you are using the external ebq 100 supply unit, connect it to an AC power outlet at this time (socket 77.1).


Fig. 76 Rear view of Leica CTR4000

- 1 AC power socket
- 2 **Microscope** socket for microscope
- 3 socket **12 V, max 100 W** for the lamp power cable of stand

Fig. 77 Rear panel of ebq 100 supply unit

- 1 AC power socket

7. Start-up

7.1 Functional Principle

Thanks to its intelligent automation, the Leica DMI5000M can be controlled using a variety of control elements.

1. Intelligent automation

- Switching between contrast methods at the touch of a button. DIC prisms, etc. are automatically positioned in the beam path.
- The microscope recognizes the selected objective and associated contrast method. The intensity (INT), aperture diaphragm (AP) and field diaphragm (FD) are always set to optimal values.
- the INT, AP and FD values are always based on the currently activated illumination axis (transmitted light or incident light).
- The INT, AP and FD values can be adjusted individually. Manual adjustments overwrite the previous settings. The current setting is stored and is retained from one session to the next when power is switched off.

2. Controls

- Rotary buttons on SmartMove (remote control) for stage (X/Y) and focus (Z) control
- Fixed function buttons on stand for INT, AP and FD, as well as for switching between transmitted-light and incident-light axis
- Variable function buttons on stand and SmartMove
These function buttons have functions suitable to the configuration of your microscope assigned to them at the factory. The functions can be reprogrammed and/or adapted to your specific requirements, however.
- Complete control of microscope and camera via software (Leica Application Suite (LAS))

7. Start-up

Note: (reset function)

The table on the following page provides an overview of the microscope functions and their controls.

The microscope can be reset to its factory default programming:

- With the stand switched off, press the top three variable function buttons on the left side of the stand.
- Switch on the power for the stand.
- Hold the buttons until the initialization is complete.
- The standard information display will now appear in the LeicaScreen (Fig. 79 and 80, p. 55).
- Switch the instrument off and back on. The settings are now saved.

7. Start-up

Function	Fixed Function buttons Stand	Variable Function buttons Stand	SmartMove Function buttons	Rotary knobs	Software
Change contrast method	-	+	+	-	+
Change transmitted light/incident-light axis	+	-	-	-	+
Change to objective	-	+	+	-	+
Teach-in parfocality	-	-	-	-	+
Change operating mode (dry/imm)	-	+	+	-	+
Illumination Manager	+	+	+	-	+
Magnification changer	+	-	-	-	+
Focusing	+	-	-	+ ¹⁾	+
Set stops	+	-	-	-	+
Go to stop	+	-	-	-	+
Change step increment (coarse/fine) -		+	+	-	+
XY stage positioning	-	-	-	+	+
Change speed	-	-	-	-	+
Stage positions (store/go to)	-	-	-	-	+
Change to reflector cube	+	(+)	+	-	+
Side and bottom port	+	(+)	+	-	+
DIC fine adjustment	+	-	-	-	+

+ always possible

(+) optional

- not possible

¹⁾ Focusing alternatively via wheels on the stand

7. Start-up

Possible assignments for variable function buttons on stand and SmartMove

For Leica DMI5000 M

Function button	Function
BF (RL)	Bright field incident light
ICR	Interference contrast, incident (reflected) light
DF (RL)	Dark field incident light
POL (RL)	Polarization incident light
CHANGE RL ●	Toggle through all incident-light contrast methods
FF	Focus Finder
BF (TL)	Bright field transmitted light
ICT	Interference contrast, transmitted light
DF (TL)	Dark field transmitted light
POL (TL)	Polarization transmitted light
PH	Phase contrast transmitted light
CHANGE TL ●	Toggle through all contrast methods
INT ↑	Increase intensity
INT ↓	Reduce intensity
AP ↑	Open aperture diaphragm (incident light)
AP ↓	Close aperture diaphragm (incident light)
FD ↑	Open field diaphragm (incident light)
FD ↓	Close field diaphragm (incident light)
SHUTTER RL	Open/close incident-light shutter
FLUO	Fluorescence (last filter cube)
CHANGE OBJ CW	Cycle through objectives clockwise
CHANGE OBJ CCW	Cycle through objectives counterclockwise
Z FINE	Activate fine focusing
Z COARSE	Activate coarse focusing
XY PRECISE	Activate precise stage
XY FAST	Activate fast stage
DRY/IMM	Switch dry/immersion
CHANGE FLT	Switch TL filter

7.2 Switching on the Microscope

- Switch on the power of the Leica electronics box at the On/Off switch (78.1). The signal lamp (78.2) is lit green when the unit is ready. All motorized microscope components will then run through an initialization phase.

Note:

If a PC is connected, switch on **the electronics box first**, and **then the computer**.

After the initialization (Fig. 79) is complete, the LeicaScreen will display the microscope's current settings (Fig. 80).

If a component has not been installed correctly, the LeicaScreen will display an error message. See Troubleshooting chapter, → p. 87.

Components such as diaphragms, condensers, light and phase rings have been pre-centered at the factory. It may be necessary to correct the centering after the microscope has been transported and assembled.

Before performing the required steps, please familiarize yourself with the LeicaScreen and the controls.

Caution!

After turning on the gas discharge lamp, the burner must be immediately adjusted. Therefore, **do not** turn on the power supply unit yet. First, work in transmitted light in order to familiarize yourself with the microscope's controls.

Fig. 78 Leica CTR4000 and CTR6000 front panel

- On/Off switch
- Signal lamp

Fig. 79
LeicaScreen
Initialization

Fig. 80
LeicaScreen
after
Initialization

IL_BF	
10x Obj.	BF
1.5x MagCh.	Σ 150x
INT 15	
AP 10 ○	FD 2 ○
80%	20%

7. Start-up

7.3 The LeicaScreen

The screen displays the microscope's current settings. The content of the display depends on the features of the individual microscope. For information on the abbreviations used, please turn to the table of abbreviations → p. 95f.

The screen has a number of areas and lines.

- Line 1: contrast method
- Line 2: objective/magnification
- Line 3: illumination/diaphragms
- Line 4: active ports
- Line 5: focus/stops

The content of the display changes according to the active function.

Pictograms


	Contrasting method
	Objective/ Magnification
	Illumination/ Diaphragms
	Ports/Eyepiece
	Focus/stops

Fig. 76 Arrangement of the function buttons – overview

- 1 Four variable function buttons
- 2 Illumination Manager
- 3 Front control panel
- 4 Focus control buttons
- 5 Three variable function buttons
- 6 SmartMove knobs
- 7 SmartMove function buttons

7.4 The Function Buttons on the Stand

A number of function buttons are located on both sides of the stand. These can be broken down into fixed and variable buttons. The variable function buttons have different functions depending on the features of the individual microscope.

Note:

Changes to the light intensity as well as aperture and field diaphragm settings are stored for the individual objectives and contrast methods.

Fixed function buttons on the left side


The **TL/IL** button (82.1) toggles between the incident-light and transmitted-light axis. The contrast method last used with a given axis is restored when switching.

The **INT** buttons (82.3) adjust the light intensity. The adjustment can be made in coarse or fine steps. Pressing both **INT** buttons at the same time toggles between coarse and fine adjustment. "Intensity fine" will appear in the display when fine adjustment is selected.

The **AP** buttons (82.2) for the aperture diaphragm and **FD** (82.4) for the field diaphragm open and close their respective diaphragms.

Fig. 82 Fixed function buttons (left side of stand)

- 1 Toggle transmitted light/incident light
- 2 Aperture diaphragm
- 3 Light intensity
- 4 Field diaphragm

7. Start-up

Variable function buttons on the stand

The variable function buttons are assigned functions at the factory that are appropriate to the features of your microscope. They are labeled accordingly. For details on button assignments, please refer to the included identification sheet.

For information on the abbreviations used, please refer to the list → p. 54.

Note:

The Leica Application Suite (LAS) software is required for changing the button assignments.

Possible functions*:

BF	SHUTTER RL
ICR	FLUO
ICT	FF
DF	CHANGE OBJ CW
POL	CHANGE OBJ CCW
PH	Z FINE
CHANGE RL ●	Z COARSE
CHANGE TL ●	XY PRECISE
INT ↑	XY FAST
INT ↓	DRY/IMM
AP ↑	CHANGE FLT
AP ↓	MEM 1-6
FD ↑	
FD ↓	

Fig. 83 Function buttons (left side of stand)

- 1 Variable function buttons
- 2 Open/close aperture diaphragm
- 3 TL/IL switching
- 4 Open/close field diaphragm
- 5 Increase/decrease light intensity



Fig. 84 Function buttons (left side of stand)

- 1 Variable function buttons

* See page 54 for abbreviations

7. Start-up

Function buttons on the front panel (Fig. 85)

 100% of the light goes to the eyepiece (85.1).

←  → Toggle function for the side ports (85.2). This function depends on the individual microscope configuration.

Note:

Switching of the top port: manual

SHUTTER Opens and closes the shutter (85.3).

 Switches between the possible magnifications of the magnification changer (85.4).

 The magnification changer is set to the magnification 1x (85.5).

CUBE

The CUBE 1 to CUBE 4 (85.6) buttons permit individual cubes to be selected directly.

Press the CUBE 2 and CUBE 3 buttons at the same time to display the assignments of the variable function buttons. To reset the display, press the buttons again or wait 3 seconds.

Focus buttons (Fig. 86)

Z↑ Moves the Z drive in the indicated direction.
Z↓

SET+Z↑ Sets the upper focus stop.

SET + Z↓ Sets the lower stop.

Fig. 85 Front control panel

- 1 100% light to eyepiece
- 2 Toggle ports
- 3 Shutter
- 4 Switch between intermediate magnifications
- 5 Intermediate magnification 1x
- 6 Selecting reflector cubes(contrast method)


Fig. 86

- 1 Focus control buttons
- 2 Open reflector drawer

7. Start-up

7.5 The SmartMove Remote Control Module

SmartMove knobs

Use the knobs 87.1 and 87.2 to move the stage in X and Y directions. The image is focused using the knob (87.3).

The height of the knobs can be adjusted to a comfortable working position by turning the wheel (87.4).

Variable function buttons on SmartMove

The variable function buttons are assigned functions at the factory that are appropriate to the features of your microscope. They are labeled accordingly. For details on button assignments, please refer to the included identification sheet.

For information on the abbreviations used, please refer to the list → p. 54.


Note:

The Leica Application Suite (LAS) software is required for changing the button assignments.

Fig. 87 SmartMove remote control module

- 1 Travel in x
- 2 Travel in y
- 3 Focus
- 4 Individual adjustment of button height
- 5 Variable function buttons (factory preset)

7.6 Incident Light

Suitable aperture and field diaphragm values have been preset for each objective. The incident light module has also been centered at the factory.

However, it may be necessary to readjust the incident light module in some cases after transporting and setting up the stand. Therefore, check the diaphragm centering.

The following procedure is provided for the incident light-bright field illumination.

- Select an objective with moderate magnification (10x-20x).
- Activate the incident-light axis with the **TL/IL** button (82.1).
- Press the **IL-BF** button to activate the bright field contrast method (one of the variable function buttons on the stand).
- Place a specimen on the stage.
- Focus the specimen using the SmartMove or the focus wheels.

- Adjust the light intensity with the **INT** buttons (82.3).


Adjusting the field diaphragm

- Close the field diaphragm with the **FD** button (82.4) or manually until the edge of the diaphragm (round or rectangular) appears in the field of view.
- If the limits of the field diaphragm are not in the center of the field of view, move the position of the field diaphragm to the center with the two centering screws (88.1) on the right side of the stand.
- Use the function buttons **FD** (82.4) to open the field diaphragm to the point that they just disappear from the field of view.
- We recommend the use of a rectangular field diaphragm when using a digital camera. Match the size of the diaphragm to the chip size of the camera.

Adjusting the aperture diaphragm

- Remove one eyepiece (e.g. right).
- Close the aperture diaphragm with the **AP** button (82.2) until the edge of the diaphragm appears in the exit pupil of the objective (aperture diaphragm plane).
- If the image is not in the center of the exit pupil, move the position of the aperture diaphragm to the center with the two centering screws (88.2).
- Open the aperture diaphragm to cover 2/3 of the field of view.

Fig. 88 Adjusting the diaphragms (incident light)
1 Adjusting screws for moving the field diaphragm
2 Adjusting screws for moving the aperture diaphragm

7. Start-up

7.7 Adjusting the Light Sources

Transmitted-light axis (TL) with lamp housing 107/2

The lamp housing 107/2 with a 12 V 100 W halogen lamp is fixed. Centering the lamp is not required.

Lamp housing 107 L for 12 V 100 W halogen lamp

The lamp can be adjusted using the screws (89.2) and the button (89.3).

- Place a sheet of white paper under the field diaphragm.
- Adjust the lamp to create an evenly bright spot on the paper.

Incident-light axis (IL) with lamp housing 106 z

- When a supply unit is used, it is turned on first.
- Activate the incident-light axis with the **TL/IL** function button. **IL** will appear on the LeicaScreen.
- Insert the lamp adjustment reflector (Fig. 90) in the filter turret in place of a filter cube. Make a note of the designation of the replaced filter cube.
- Turn the reflector into the beam path. The reflector is correctly positioned when the LeicaScreen shows the designation of the replaced filter cube.

Fig. 89 Lamp housing 107 L

- 1 Locking screws for housing
- 2 Screw for vertical adjustment
- 3 Button for horizontal adjustment
- 4 Collector focusing



Fig. 90 Reflector cube for lamp adjustment

**Caution!**

Never look directly into the beam path!
Beware of the glare hazard when switching
to reflector BF or Smith!

**Caution!**

Light sources pose a potential irradiation
risk (glare, UV-radiation, IR-radiation).

In the lamp housing 106z, the direct image of the filament (in halogen lamps) or the arc (in gas discharge lamps) and its reflection are focused separately and adjusted in relation to one another.

An adjustment window is located on the right side of the microscope (3.4, p. 17) in which the light source is visible.


Adjust the lamp as follows while observing the light source in the adjustment window.

Centering the Hg 100 W and Xe 75 W mercury lamps

- The adjustment window shows the direct image of the arc and its mirror image. These are generally not in alignment with one another.
- Focus the direct image with the collector (91.6).
- Use the adjusting buttons to pivot the arc's mirror image on the rear side of the lamp housing (91.2, 91.4,) to the side or completely out of the beam path. The arc's focused image remains visible (Fig. 92).
- Use the adjusting buttons (91.1 and 91.5) to place the direct arc image in the middle of the centering plane, whereby the bright tip of the arc, the focal spot, should lie slightly outside the center (Fig. 93).

Fig. 91 Lamp housing 106z L

- 1 Lamp adjustment, vertical
- 2 Vertical reflector adjustment
- 3 Focusing the reflector image
- 4 Horizontal reflector adjustment
- 5 Lamp adjustment, horizontal
- 6 Collector focusing
- 7 Screw

7. Start-up

- Then pivot the arc's mirror image with the adjusting knobs (91.2) and (91.4) and focus it using the reflector (91.3).
- Use the adjusting knobs (91.2) and (91.4) to orient the mirror image symmetrically to the direct image (Fig. 94).
The V-shaped irradiation of the direct image and mirror image arcs can be superimposed.

Caution!

The bright tips of the arcs, the focal spots, must never be projected onto each other, as this results in a danger of explosion by overheating.

Caution!

The structure of the arc can no longer be made out clearly in lamps that have been in service for a long time. The image and mirror image can no longer be superimposed exactly. In this case, align both images.

- Using the collector, defocus the image with the knob (91.6) until the arc image and mirror image are no longer recognizable and the image is homogeneously illuminated.
- Replace the lamp adjustment reflector with the original filter cube.

Fig. 92 Direct arc image focused but not centered (in reality, the image is less focused)


Fig. 93 Direct arc image in target position (in reality, the image is less focused)



Fig. 94 Direct arc image and mirror image in target position (in reality, the image is less focused)

8. Operation

8.1 Switching on

Switch on the power of the Leica electronics box at the On/Off switch (95.1). The signal lamp (95.2) is lit green when the unit is ready.

Note:

If you have connected a PC, switch the **electronics box on first, then open the Leica Application Suite (LAS)**.

All motorized microscope components will then run through an initialization phase.

Note:

In the case of faulty initialization ("Init Error" message on LeicaScreen), see Troubleshooting chapter, → p. 87ff.

All of the user's previous settings are restored during the initialization.

Caution:

The focal position and the lower stop are also retained from one session to the next when power is switched off.

After the initialization is complete, the LeicaScreen will display the status screen with the microscope's current settings. Fig. 97 is an example.

Fig. 95
Front side
Leica CTR4000 and CTR6000
1 On/Off switch
2 Signal lamp

Fig. 96
LeicaScreen
Initialization

Fig. 97
LeicaScreen
after
Initialization

	IL_BF	
	10x Obj.	BF
	1.5x MagCh.	Σ 150x
	INT 15	
	AP 10 ○	FD 2 ○
	80%	20%

8. Operation

8.2. Contrast Methods

All the contrast methods of the Leica DMI5000 M can be selected and controlled via the variable function buttons and the Leica Application Suite (LAS). The only exceptions are methods that involve components requiring manual control (e.g. systems with manual analyzers). The following section describes the use of the function buttons on the stand. For instructions on the use of the software, please refer to the separate manual.

8.2.1 Bright field (RL/TL)

- Use the **TL/IL** function button to switch to incident light (TL).
- Select the **BF** (bright field) contrast method. Press the variable button to which you have assigned the **BF** function. Alternatively: press the variable button **CHANGE RL** – this will toggle through all of the contrast methods supported by the current configuration. (For details on button assignments, please see the identification sheet.) **BF** will appear on the LeicaScreen.
- Insert an incident light specimen.
- Bring an appropriate objective into place.
- Focus the image with the knob on the SmartMove or the focusing wheel and adjust the intensity with the **INT** function buttons.

Fig. 98 Function buttons (left side of stand)

- 1 Variable function buttons
- 2 Open/close aperture diaphragm
- 3 TL/IL switching
- 4 Open/close field diaphragm
- 5 Increase/decrease light intensity



Fig. 99 Function buttons (right side of stand)

- 1 Variable function buttons

8. Operation

The procedure is similar for bright-field illumination in transmitted light.

- Use the **TL/IL** function button to switch to transmitted light (TL).
- Select the **BF** (bright field) contrast method
Press the variable button to which you have assigned the **BF** function.
Alternatively: press the variable button **CHANGE TL**  – this will toggle through all of the contrast methods supported by the current configuration.
(For details on button assignments, please see the identification sheet.)
BF will appear on the LeicaScreen.
The motorized condenser will now move to the bright field position.
The bright-field position of the 4-position reflector turret disk will now be selected.

8.2.2 Dark field (RL/TL)

- Use the **TL/IL** function button to switch to incident light (RL).
- Select the **DF** (dark field) contrast method by pressing the variable button **DF**.
Alternatively: press the variable button **CHANGE RL**  .
(For details on button assignments, please see the identification sheet.)
DF will appear on the LeicaScreen.
The dark-field position of the 4-position reflector turret disk will now be selected.
- Insert a incident light specimen.
- Bring an appropriate objective into place.

- Focus the image with the knob on the SmartMove or the focusing wheel and adjust the intensity with the **INT** function buttons.

The procedure is similar for bright-field illumination in transmitted light.

- Use the **TL/IL** function button to switch to transmitted light (TL).
- Select the **DF** (bright field) contrast method
Press the variable button to which you have assigned the **DF** function.
Alternatively: press the variable button **CHANGE TL**  – this will toggle through all of the contrast methods supported by the current configuration.
(For details on button assignments, please see the identification sheet.)
DF will appear on the LeicaScreen.
The motorized condenser will now switch to the dark field ring.

Note:

When selecting the dark field method, the aperture diaphragm is opened fully.

8. Operation

8.2.3 Polarization (RL/TL)

- Use the **TL/IL** function button to switch to incident light (IL) or transmitted light (TL).
- Select the **POL** (polarization) contrast method by pressing the variable button **POL**.
Alternatively: press the variable buttons **CHANGE RL** or **TL**.
(For details on button assignments, please see the identification sheet.)
POL will appear on the LeicaScreen.

Manual method:

- Incident light: move the polarizer into the beam path.
Transmitted light: move the polarizer on the condenser into the beam path.
- Insert the analyzer into the right side of the stand until it clicks into position (Fig. 100).
- Bring the polarizer and analyzer into cross position until they reach maximum darkness.
- Place a specimen on the stage and select a suitable objective.

Motorized method:

- If the microscope is equipped with the relevant components, the **POL** position of the 4-position reflector turret disk will be selected automatically when the **POL** contrast method is selected with incident illumination. The Pol or ICR reflector is positioned in the beam path. With transmitted illumination, the polarizer is automatically positioned in the condenser and the analyzer cube automatically positioned in the beam path.

Combined methods:

- The Leica DMI5000 M microscope supports the combined use of automatic and purely manual components.

Fig. 100 Inserting the analyzer

8.2.4 Differential Interference Contrast (RL/TL)

- Use the **TL/IL** function button to switch to incident light (IL) or transmitted light (TL).
- Select the **DIC** contrast method
Press the **ICR** or **ICT** variable buttons.
Alternatively: press the variable buttons **CHANGE RL** or **TL**.
(For details on button assignments, please see the identification sheet.)
DIC will appear on the LeicaScreen.
- Incident light:
The ICR reflector cube on the 4-position reflector turret disk and the suitable objective prism are automatically positioned in the beam path. The corresponding objective prism and the analyzer cube are also positioned automatically.
Transmitted light:
The polarizer in the condenser and the suitable condenser prism are automatically positioned in the beam path. The corresponding objective prism and the analyzer cube are also positioned automatically.
- Place a DIC specimen on the stage.
- Bring an appropriate objective into place.
- Focus the image with the knob on the SmartMove or the focusing wheel and adjust the intensity with the **INT** function buttons.
- Use the knurled wheel below the objective turret for fine adjustment (Fig. 101).

Manual alternative:

- Incident light: move the polarizer into the beam path manually.
Transmitted light: move the polarizer on the condenser into the beam path manually.
- Insert the analyzer manually into the right side of the stand until it clicks into position (Fig. 100).
- Use the knurled wheel below the objective turret for fine adjustment (Fig. 101).

Fig. 101 DIC disk with knurled wheel for fine adjustment

8. Operation

8.2.5 Fluorescence

- Use the **TL/IL** function button to switch to fluorescence **FLUO**.
- Place a specimen on the stage and select a suitable objective.
- The current fluorescence filter cube will be displayed on the LeicaScreen.
- You may protect your specimen from fading by closing the incident light shutter.
To do so, press the **SHUTTER** button (85.3, p. 59) on the front panel.
The following pictogram will appear on the LeicaScreen:

► Changing the fluorescence filter cube

- Function buttons on the front panel: **CUBE 1** and **CUBE 4**
- Leica Application Suite (LAS) Software

- Focus the image with the knob on the SmartMove or the focusing wheel and adjust the intensity with the **INT** function buttons.

8.2.6 Phase Contrast (TL)

- Use the **TL/IL** function button to switch to transmitted light (TL).
- Select the **PH** (phase contrast) contrast method by pressing the variable button **PH**.
Alternatively: press the variable button **CHANGE TL** .
(For details on button assignments, please see the identification sheet.)
PH will appear on the LeicaScreen.
The motorized condenser will now switch to the correct light ring.
- Insert a transmitted light specimen.
- Rotate an appropriate objective into place.
Objectives that are suitable for phase contrast are engraved with PH.
- Focus the image with the knob on the SmartMove or the focusing wheel and adjust the intensity with the **INT** function buttons.

Note:

When selecting the phase contrast method, the aperture diaphragm is opened fully and can not be adjusted.

8.3 Focusing

Note:

In **DMI5000 M instruments with motorized focus**, the parfocality teach-in has already been performed at the factory. However, it may be necessary to perform another teach-in after installing the objectives when setting the microscope up. We recommend checking parfocality before setting the stops and performing a teach-in with the Leica Application Suite (LAS).

Focusing the image

The image is focused using the knob (87.3 p. 60) on the SmartMove remote control module

Alternatively, use the focus wheels on either side of the stand.

The current Z position is shown on the LeicaScreen. In the case of motorized stages, the Z drive will travel to its lowest position prior to the stage initialization when switching the microscope on.

The focus buttons **Z↑** and **Z↓** on the right side of the stand (Fig. 102) permit fast focusing or lowering of the objectives.

Fig. 102

1 Focus control buttons

Setting Low Threshold

Set the lower focus threshold by pressing and holding the **SET** button and pressing the **Z↓** button as well.

The display will show ▼ .

Pressing the button combination again will delete the low threshold.

The display will show ▼ .

The lower threshold can also be set using the Leica Application Suite (LAS).

The **lower threshold** is the same for all objectives and can not be traversed.

In addition, the **focus position** can also be set.

To do so, press and hold the **SET** button and press the **Z↑** button as well.

The display will show ✕ .

Pressing the button combination again will delete the focus position.

The display will show ✕ .

The focus position can also be set using the Leica Application Suite (LAS).

Set the focus position for the dry objective at the highest magnification. The focus positions will then be set automatically for all other objectives, taking parfocality and working distances into account.

► Set the stops via

- fixed function buttons on stand
- Leica Application Suite (LAS) Software

8. Operation

Summary of pictograms

- ▼ lower threshold not set
- ▼ lower threshold set
- ✕ focus position not set
- ✕ focus position set

Going to the Lower Threshold

Go to the lower threshold by pressing and holding the **Z↓** button.

Go to the focus position by pressing and holding the **Z↑** button.

These functions can be controlled via software.

► Go to stops via

- fixed function buttons on stand
- Leica Application Suite (LAS) Software

Note:

When going to the stops with the **Z↑** and **Z↓** buttons, hold the button until the stop has been reached.

Setting the step increments

It is possible to toggle between **Fine** and **Coarse** step increments.

The **Fine** value varies to suit the current objective. Suitable values have been predefined. The assignments can be changed with the Leica Application Suite (LAS).

When selecting **Coarse**, the positioning speed is the same for all objectives. **Coarse** corresponds to the maximum speed.

Note:

The assignment of a specific step increment to an objective not only applies to the Z drive, but also to the step increments assigned to the stage when **Precise** (→ p. 84) is selected.

► Switch between Fine and Coarse via

- variable function buttons on stand and SmartMove
- Leica Application Suite (LAS) Software

8.4 Tubes

Note:

Close any unused tube openings, as otherwise stray light can interfere with observation.

Adjusting the viewing distance

- Adjust the viewing distance of the eyepieces so that a congruent total image is seen (Fig. 103).

Adjusting the viewing angle

- Ergotubes feature a tilting binocular section for a 30°–45° viewing angle adjustment range.

Beam splitting in photo tubes

The beam splitting is set manually by pulling out a control bar.

Button		Observation	Photo
VIS		100 %	0 %
50/50		0 %	100 %

► Light distribution via

- manual control bar

Port selection

The

button on the front control panel switches 100% of the light to the eyepieces.

Use the

button, also on the front control panel, to select the side ports.

Depending on the configuration, the screen will now display

- the active port (right or left) and
- the percentage of light going to the port (100%, 80%, 50%).

► Select ports via

- fixed function buttons on stand (side ports)
- Leica Application Suite (LAS) Software
- manual action (top port)

Fig. 103 Tube setting

8. Operation

8.5 Optical Equipment

Objectives that have been optimized for incident-light specimens without coverslips (cover-slip thickness 0) and incident-light methods are preferable for metallographic examinations.

These objectives are also fully usable for trans-mitted-light microscopy if the transmitted-light specimens have not been coverslipped. Objectives with low magnification and a numeri-cal aperture of less than 0.25 can also be used with coverslipped transmitted-light specimens without losses.

Please observe the objective identification:

Objective identification (examples):

N-PLAN	10 x	0.25		×	-*	A
HC PL FLUOTAR	50 x	0.80	BD	×	0	D
Correction class	Image scale	Numerical aperture	Bright and dark field objective	Infinite Optics	Coverslip correction	ICR prism type

Eyepiece identification:

10x = magnification
/20 = field of view (mm)
M = adjustable eyelens
 = for eyeglass wearers (antiglare hood removable or folding eyecup)

* - = designed for specimens with or without coverslips
0 * - = designed for specimens without coverslips

The DMI5000 M microscope is based on the tube length \times (infinite) and a tube lens reference focal length of $f = 200$ mm.

Fig. 104 Optical equipment






- 1 HC PL FLUOTAR BD objective series
M32 x/0.75 thread
- 2 Eyepieces HC PLAN 10 x/20  and  M
Eyepieces HC PLAN 10 x/22  and  M
Eyepieces HC PLAN 10 x/25  and  M

Fig. 105 Optical equipment

- 1 N-PLAN objective series
for bright field, M25 x 0.75 thread
- 2 Eyepieces HC PLAN 10 x/20  and  M
Eyepieces HC PLAN 10 x/22  and  M

8. Operation

In the DMI5000 M Leica objectives with infinite optics must be used.

Note:

Older Leica objectives with finite optics correction may only be used with bright-field versions with 32/RMS spacer rings.

Depending on the year of manufacture of the objective, the engraved magnification values may vary by a factor of 0.8x.

Objectives for bright-field incident light, polarization contrast, interference contrast*

Note: Bright-field objectives with M25 x 0.75 mm thread or RMS thread need spacer rings for the M32 x 0.75 mm thread of the objective turret (Fig. 24, p. 30)

with N PLAN objectives with M25 x 0.75 thread

N PLAN	2.5 x/0.07	× /-	FAA* 11.20 mm	Order No. 506083
N PLAN	5 x/0.12	× /-/A	FAA* 14.00 mm	Order No. 506087
N PLAN	10 x/0.25	× /-/A	FAA* 5.80 mm	Order No. 506084
N PLAN	20 x/0.40	× /0/D	FAA* 1.10 mm	Order No. 566026
N PLAN	50 x/0.75	× /0/D	FAA* 0.37 mm	Order No. 566027
N PLAN	100 x/0.90	× /0/D	FAA* 0.27 mm	Order No. 566028
Spacer ring	32/25 6 x			Order No. 561003

with PL FLUOTAR objectives with M25 x 0.75 thread

PL FLUOTAR	1.6 x/0.05	× /-	FAA* 1.54 mm	Order No. 566010
PL FLUOTAR	2.5 x/0.07	× /-	FAA* 9.20 mm	Order No. 567010
HC PL FLUOTAR	5 x/0.15	× /-/D	FAA* 12.00 mm	Order No. 506504
HC PL FLUOTAR	10 x/0.30	× /-/D	FAA* 11.00 mm	Order No. 506505
HC PL FLUOTAR	20 x/0.50	× /0/D	FAA* 1.27 mm	Order No. 566500
HC PL FLUOTAR	50 x/0.80	× /0/D	FAA* 0.50 mm	Order No. 566501
HC PL FLUOTAR	100 x/0.90	× /0/D	FAA* 0.30 mm	Order No. 566502
Spacer ring	32/25 7 x			Order No. 561003

High-aperture oil-immersion objective for high resolution

PL FLUOTAR	100 x/1.30	× / Oil D	FAA* 0.13 mm	Order No. 506046
Immersion oil,	DIN 58884	10 ml		Order No. 513787
Spacer ring	32/25 x			Order No. 561003

* FAA = free working distance

8. Operation

PL Apo series with RMS thread

PL APO	50 x/0.90	× /0/C	FAA* 0.28 mm	Order No. 567034
PL APO	100 x/0.95	× /0/C	FAA* 0.16 mm	Order No. 567023
PL APO	150 x/0.95	× /0/C	FAA* 0.20 mm	Order No. 567042
PL APO	250 x/0.95	× /0/	FAA* 0.24 mm	Order No. 767001
Spacer ring	32 x			Order No. 562281

Objectives with especially long working distances (RMS)

PL FLUOTAR L	50 x/0.55	× /0/C	FAA* 8.00 mm	Order No. 767002
PL FLUOTAR L	100 x/0.75	× /0/	FAA* 4.70 mm	Order No. 767000
Spacer ring	32 x			Order No. 562281

Objectives with long working distances with 1.8 mm quartz glass coverslip correction (M25 x 0.75)

PLAN H	20 x/0.40	× /1,8Q B	FAA* 12.60 mm	Order No. 566003
PLAN H	40 x/0.60	× /1,8Q B	FAA* 7.10 mm	Order No. 566004
Spacer ring	32/25 x			Order No. 561003

Objectives for bright-field and dark-field incident light, polarization contrast, interference contrast* with M32 x 0.75 mm thread*

Objectives for bright-field and dark-field with M32 x 0.75 mm thread screw straight into the turret without spacer rings.

N PLAN Objectives

N PLAN	5 x/0.12 BD	× /-/A	FAA* 13.20 mm	Order No. 566016
N PLAN	10 x/0.25 BD	× /-/A	FAA* 5.20 mm	Order No. 566005
N PLAN	20 x/0.40 BD	× /0/D	FAA* 1.10 mm	Order No. 566029
N PLAN	50 x/0.75 BD	× /0/D	FAA* 0.37 mm	Order No. 566030
N PLAN	100 x/0.90 BD	× /0/D	FAA* 0.30 mm	Order No. 566031

PL FLUOTAR Objectives

HC PL FLUOTAR	5 x/0.15 BD	× /-/D	FAA* 12.20 mm	Order No. 566506
HC PL FLUOTAR	10 x/0.30 BD	× /-/D	FAA* 11.00 mm	Order No. 506503
HC PL FLUOTAR	20 x/0.50 BD	× /0/D	FAA* 1.27 mm	Order No. 566507
HC PL FLUOTAR	50 x/0.80 BD	× /0/D	FAA* 0.50 mm	Order No. 566504
HC PL FLUOTAR	100 x/0.90 BD	× /0/D	FAA* 0.30 mm	Order No. 566505

PL Apo Objectives

PL APO	50 x/0.90 BD	× /0/C	FAA* 0.34 mm	Order No. 567013
PL APO	100 x/0.95 BD	× /0/C	FAA* 0.26 mm	Order No. 567014
PL APO	150 x/0.95 BD	× /0/C	FAA* 0.25 mm	Order No. 567015

Objectives with long working distance

PL FLUOTAR L	20 x/0.40 BD	× /0/C	FAA* 11.10 mm	Order No. 766001
PL FLUOTAR L	50 x/0.55 BD	× /0/C	FAA* 8.10 mm	Order No. 766000

* FAA = free working distance

8. Operation

The following eyepieces are recommended to obtain the standard magnification:

Eyepiece HC PLAN S 10 x/25 M Order No. 507808

Eyepiece HC PLAN 10 x/22 M 2 x Order No. 507804

Alternative eyepieces for the above objectives:

Eyepiece HC PLAN 10 x/20 Order No. 507801

Eyepiece HC PLAN 10 x/20 M Order No. 507802

Other eyepieces outside the standard magnification:

Eyepiece HC PLAN 12.5 x/16  M Order No. 506515

Eyepiece 16 x/14 B, adjustable, Order No. 445301 Leica Heerbrugg

Eyepiece 25 x/9.5 B, adjustable, Order No. 445302 Leica Heerbrugg

Non-HC series eyepieces from Leica Heerbrugg require a spacer ring, order No. 506808

8. Operation

8.6 Graticules for Eyepieces

Graticules for length measurement, grain and particle size determination

To calibrate the graticules, we recommend:

Incident light object micrometer
1 mm = 100 divisions Order No. 563011

The following graticules are available:

- Graticule
10 mm/100 divisions Order No. 506950
- Graticule 10 mm/100 divisions
with crosshair Order No. 506952
- Graticule for standard series
and Snyder-Graf method Order No. 566950
- Graticule ASTM-E-112
Grain size determination Order No. 566951
- Graticule with grid
10 x 10 x 0.1 mm Order No. 506954
- Graticule with grid
10 x 10 x 1 mm Order No. 506955
- Graticule with
crosshair Order No. 506953
- Graticule
10 mm/200 divisions Order No. 506951
- Format graticule F6
for microphoto (for MPS with
photo eyepiece HC 10x) Order No. 506951
- Format graticule F7
for microphoto (for DMLD with
photo eyepiece HC 10x) Order No. 506962
- Format graticule F8
for microphoto
(for DMLD and MPS with
photo eyepiece HC 12.5x) Order No. 506963

Graticule 10 mm/100 T (Fig. 106)

When using gratitudes with a scale of 10 mm/100 divisions, one division corresponds to approximately the following lengths in the specimen plane:

0.02 mm =	20 μm	with	5x objective
0.01 mm =	10 μm	with	10x objective
0.005 mm =	5 μm	with	20x objective
0.002 mm =	2 μm	with	50x objective
0.001 mm =	1 μm	with	100x objective

Graticule for standard series/Snyder Graff (Fig. 107)

The graticule for standard series and the Snyder Graff method contains a central circle that marks an image detail of 80 μm at standard magnification.

The 80 μm dia. image detail corresponds to the grain images of standard series charts and permits the exact comparison of grain sizes.

The central circle can also be used for determining grain sizes with the circle segment or area count methods.

The graticule also features a baseline for the Snyder Graff method.




This method involves counting the number of grains intersected by the reference line and then calculating the average grain size from several measurements.

ASTM-E-112 graticule (Fig. 108)

The ASTM-E-112 is divided up into 8 segments with labeled grain size pictures.

The pictures conform to grain size plates No. 1A and 1B of the ASTM-E-112 standard.

For reasons of space, the graticule contains 8 representative circular grain size images from a total of 22 of the above-mentioned plates 1A and 1B.

Fig. 106 Graticule 10 mm/100 T**Fig. 107** Graticule for standard series and Snyder Graff**Fig. 108** Graticule ASTM-E-112

8. Operation

Inserting the graticules:

The graticules have a diameter of 26 mm and can only be inserted in the eyepieces HC PLAN 10 x/20/22/25 with type M adjustable eyelens. The second eyepiece in the tube should also have a type M adjustable eyelens.

Caution:

Make sure there are no dust particles or fingerprints on the graticule when inserting it, as these will show up on the microscope image.

HC PLAN 10 x/20 M and HC PLAN 12.5 x/16 M eyepieces

Unscrew the retainer ring on the underside of the eyepiece (109.8a), using a rubber cloth if necessary.

HC PLAN 10 x/22 and HC PLAN 10 x/25 eyepieces

Unscrew the lower part of the eyepiece (109.7). Unscrew the slitted retainer ring (109.8b) out of the lower part using a blunt blade. Insert the graticule with the coated side pointing downward toward the objective. Lettering must be shown the right way around. Screw the retainer ring and lower part of the eyepiece back in.

Fig. 109 Eyepieces

- 1–4** Eyepieces for users without eyeglasses (antiglare eyecup 10 in place or extended)
- 5** FOTO eyepiece
- 6** 10 x/25M eyepiece disassembled - upper part
- 7** lower part, unscrewed (also applies to 10 x/22M, 12.5 x/16M, but not to 10 x/20 and 10 x/20M)
- 8a,b** Retaining ring for eyepiece graticules, removable
- 9** Eyepiece graticule
- 10** Antiglare hood, removable for observation with eyeglasses (foldable on 10 x/20 and 10 x/22 eyepieces, insertion and removal item 8a and 8b). The design of the 12.5 x/16M eyepiece essentially corresponds to that of the 10 x/25 M eyepiece.

8.7 Eyepieces

Note:

The eyepiece's aperture protector must be removed or folded back during microscopy while wearing eyeglasses.

We recommend removing bifocals and spectacles with progressive-addition lenses when using the microscope.

- For the adjustable tubes with documentation output, choose the 100% VIS position.

Eyepieces with inlaid reticle

- Focus the reticle by adjusting the eyelens.
- Focus on the object through this eyepiece.
- Then, close that eye and focus on the object by adjusting only the second ocular.

Correction for vision problems

- With your right eye, look through the right eyepiece and bring the specimen into sharp focus.
- Then, with your left eye, view the same specimen and rotate the left eyepiece tube until the object is brought into sharp focus. Do not change the Z position in the process!

Note:

We recommend running a teach-in via the Leica Application Suite (LAS) software when using eyepieces not included in the scope of delivery. This will ensure that the total magnification shown in the LeicaScreen is correct.

8.8 Objectives

Changing objectives

The objectives can be selected with the function buttons on the stand or the SmartMove, or by manually turning the objective turret. When changing objectives manually, please ensure that the turret clicks into position.

The positions of the objectives in the objective turret have been specified at the factory and must be observed when installing the objectives.
(also see Objectives →, p. 30).

When selecting an objective, the microscope automatically selects:

- the optimal setting for the field diaphragm
- the optimal setting for the aperture diaphragm
- the light intensity for the current contrast method

The objective magnification and total magnification are displayed on the LeicaScreen.

8. Operation

- For **immersion objectives** use the appropriate immersion medium.

OIL: only use optical immersion oil according to DIN/ISO standards.
Cleaning → p. 91

W: Water immersion.

IMM: Universal objective for water, glycerol, oil immersion.

Caution!

Follow safety instructions for immersion oil!

► Select objectives via

- variable function buttons on stand and SmartMove
- Leica Application Suite (LAS) Software
- Manual selection

Changing the operating modes "dry" (DRY) and "immersion" (IMM)

(for Leica DMI5000 M with motorized focus)

Each objective is assigned to a specific objective category:

- 1) Dry objectives (DRY)
- 2) Immersion objectives (IMM)

Note:

It is possible to use objectives for both operating modes.

The mode can be assigned in the Leica Application Suite (LAS).

Changing the operating mode

- First, select the operating mode (Imm or Dry) using the function buttons.
The operating mode may also be selected in the Leica Application Suite (LAS).
- The objective turret is lowered to its bottom stop. This is to permit the application of the immersion liquid when changing from a dry to an immersion objective. It also permits the removal of the liquid when changing to dry mode.
The current objective remains in the beam path.
- Next, press the button for the objective you intend to use.

Note:

If the **Imm** or **Dry** operating mode buttons are pressed accidentally, the original mode can be restored by pressing the appropriate button.

► Change operating mode via

- variable function buttons on stand and SmartMove
- Leica Application Suite (LAS) Software

8. Operation

Note:

When replacing objectives, you must perform a teach-in for the new objectives in the Leica Application Suite (LAS). A parfocality teach-in must then be performed.

Note:

For lockable immersion objectives lock these by pushing the front part upwards until it stops (approx. 2 mm). Then, after a gentle turning motion to the right, the objective is locked.

For objectives with correction mountings turn the knurl to adjust the objective to the thickness of the cover glass.

Color coding of objectives

The magnification of each objective is indicated by a color ring in accordance with DIN/ISO standards:

100x 125x 150x 160x	63x	40x 50x	25x 32x	16x 20x	10x	6.3x	4x 5x	2.5x	1.6x
white	dark-blue	light-blue	dark-green	light-green	yellow	orange	red	brown	gray

Immersion objectives are marked by an additional, lower color ring.

black oil or Imm (universal objective for oil, water or glycerin)

white water

orange glycerin

The various engraved markings of the objectives provide information on their applications:

black or dark blue bright field objectives, strain-free
green phase contrast objectives, strain-free

8. Operation

8.9 Stages and Specimen Displacement

Specimen displacement using SmartMove

The positioning of the stage is controlled by the knobs (110.1, 110.2) on the SmartMove remote control module.

Setting the step increments

The positioning speed of the stage can be varied by switching between the **Fast** and **Precise** step increments.

When selecting **Fast**, the positioning speed is the same for all objectives.

The **Precise** speed varies to suit the current objective.

- ▶ **Switch between Precise and Fast via**
 - ▶ variable function buttons on stand and SmartMove
 - ▶ Leica Application Suite (LAS) Software

Storing and restoring stage positions


A variety of stage positions can be stored temporarily in the Leica Application Suite (LAS). The XY position is stored, not the Z position.

In addition to a loading position (Load), 5 stage positions can be set temporarily. When switching the microscope on, the stage will travel to a previously-defined starting position.

- ▶ **Temporarily store and restore stage positions via**
 - ▶ Leica Application Suite (LAS) Software

Fig. 110 SmartMove remote control module

- 1 Travel in x
- 2 Travel in y
- 3 Focus (z)
- 4 Individual adjustment of button height
- 5 Variable function buttons (factory preset)

8.10 Magnification Changer

A motorized magnification changer can be used optionally. The following magnification factors are available (3 factors maximum):

1x ; 1,5x ; 1,6x; 2x

The selected factor is displayed on the LeicaScreen and in the relevant field of the Leica Application Suite (LAS), and is taken into account when calculating the total magnification.

Pressing the left button (111.1) switches between the possible magnification factors; pressing the right button selects the factor 1x.

► Change magnification via

- fixed function buttons on stand
- Leica Application Suite (LAS) Software

8.11 Light Sources

- Adjust the intensity with the function buttons (112.4). The **INT** function buttons are always assigned to the currently active transmitted light (TL) or incident light (IL) axis.

- For TL and IL:

The setting can be made in coarse and fine steps. Pressing both **INT** (112.4) buttons at the same time toggles between coarse and fine adjustment. The light intensity displayed on the LeicaScreen changes accordingly.

Coarse adjustment: 0–20

Fine adjustment: 0–255


- The intensity is individually adjusted and stored for each objective and contrast method.

► Adjust intensity via

- fixed function buttons on stand
- variable function buttons on stand and SmartMove
- Leica Application Suite (LAS) Software

Fig. 111 Front control panel

1 Function buttons for magnification changer

8. Operation

8.12 Aperture and Field Diaphragm

Both diaphragms have been set to suitable values for the current objective and contrast method at the factory.

- The diaphragms can be adjusted at any time with the **AP** (aperture diaphragm) (112.1) and **FD** (field diaphragm) (112.3) function buttons. The values displayed on the LeicaScreen change accordingly. The function buttons are assigned to the currently active transmitted light (TL) or incident light (IL) axis.

► Adjust diaphragms via

- fixed function buttons on stand
- variable function buttons on stand and SmartMove
- Leica Application Suite (LAS) Software

Caution:

The old values will be overwritten by the current ones!

Caution:

The aperture diaphragm is fully opened when using **DF**.

Note:


Focus Finder (FF)

The Focus Finder is designed to focus when working at high magnifications, with reflective surfaces, etc.

When pressing the function button to enable the Focus Finder, a smaller field diaphragm appears. Focus the field diaphragm. By focusing the diaphragm, you are also setting the focal plane of the specimen surface.

Fig. 112 Fixed function buttons, left side of stand

- 1 Aperture diaphragm
- 2 Transmitted light/incident light
- 3 Field diaphragm
- 4 Light intensity

9. Troubleshooting

Problem	Cause/Remedy
Stand	
The microscope does not respond upon being switched on.	<ul style="list-style-type: none"> ▶ Ensure that the AC outlet has power. ▶ Ensure that the electronics box is connected to an AC outlet. ▶ Check the cable connections. ▶ Inform Service and have the supply unit fuse checked.
Illumination	
The image is completely dark.	<ul style="list-style-type: none"> ▶ Open the shutter (→ p. 59). ▶ Check the connections of the lamp housings on the microscope (transmitted/incident light) ▶ Ensure that the lamps are connected to the power supply and are not defective. ▶ Inform Service and have the ebq 100 supply unit fuse checked.
The image is unevenly or not uniformly illuminated.	<ul style="list-style-type: none"> ▶ Remove all unneeded filters from the light path. ▶ Center the lamp (→ p. 62ff). ▶ Replace the old lamp (→ p. 23ff, 42f).
The illumination flickers.	<ul style="list-style-type: none"> ▶ Be sure that there is no loose connection at the power supply. ▶ Replace the old lamp (→ p. 23ff, 42f).
The fluorescence lamp does not illuminate immediately upon being switched on.	<ul style="list-style-type: none"> ▶ Hot Hg lamps should cool down before switching on again. ▶ The ebq 100 must be switched-on repeatedly. ▶ Replace the old lamp (→ p. 25ff).

9. Troubleshooting

Problem	Cause/Remedy
Bright field	
The specimen can not be brought into focus.	<ul style="list-style-type: none">▶ Use the correct immersion medium.▶ Ensure that you are using an objective without coverslip correction.
Dark Field	
No definite DF contrast is possible.	<ul style="list-style-type: none">▶ Be sure that a DF objective is being used.▶ Light intensity too low.▶ Open aperture diaphragm.
The image is unevenly or not uniformly illuminated.	<ul style="list-style-type: none">▶ The magnification is too weak. Use a higher magnification.▶ Remove the condenser head or condenser lenses.▶ Open aperture diaphragm.
Undesirable stray light.	<ul style="list-style-type: none">▶ Clean the specimen and neighboring lenses (→ p. 91).
Phase contrast	
No phase contrast is possible.	<ul style="list-style-type: none">▶ The specimen is too thick.▶ Inform the service department.

9. Troubleshooting

Problem	Cause/Remedy
Polarization	
No polarization contrast is possible.	<ul style="list-style-type: none">▶ Bring the polarizer and analyzer into cross position until they reach maximum darkness (without specimen). (→ p. 33ff, 68).
Differential Interference Contrast	
No differential interference contrast is possible	<ul style="list-style-type: none">▶ The specimen is too thick or too thin.▶ Embedding medium or specimen are of birefringent material. Rotate the specimen.▶ The difference in the refractive indices of the specimen and the embedding medium is too small.▶ The cover glass is too thick.▶ Check the Koehler illumination.▶ Bring the polarizer and analyzer into cross position until they reach maximum darkness (without specimen).▶ Check whether the suitable condenser prism and corresponding objective prism are selected (manual alternative → p. 69).▶ Ensure that the IC prisms are correctly seated (→ p. 28f).
Fluorescence	
The image is completely dark (no fluorescence).	<ul style="list-style-type: none">▶ Open the shutter (→ p. 59).▶ Select the incident-light axis (IL) (→ p. 57).▶ Insert a new lamp (→ p. 25ff).▶ Check whether the correct path has been selected at the mirror housing.
The fluorescence is too weak.	<ul style="list-style-type: none">▶ Center the lamp (→ p. 62ff).▶ Insert a new lamp (→ p. 25ff).

9. Troubleshooting

Problem	Cause/Remedy
---------	--------------

LeicaScreen

Init Error!

- ▶ Check the cable connections.
- ▶ Check whether the cover of the filter disk has clicked into position.
- ▶ Check the installed objectives, filter cubes, etc.
- ▶ Switch the microscope off and back on.
- ▶ Inform the service department.

10. Care of the microscope

Caution!

Unplug the power supply before performing cleaning and maintenance work! Protect electrical components from moisture!

Microscopes in warm and warm-damp climatic zones require special care in order to prevent the build up of fungus.

The microscope should be cleaned after each use, and the microscope optics should be kept strictly clean.

10.1 Dust Cover

Note:

To protect against dust, cover the microscope and accessories with the dust cover after each use.

Caution!

Let lamps cool down before covering the stand with a dust cover. The dust cover is not heat-resistant. In addition condensation water may occur.

10.2 Cleaning

Caution:

Residual fiber and dust can create unwanted background fluorescence.

Cleaning coated parts

Dust and loose dirt particles can be removed with a soft brush or lint-free cotton cloth.

Stubborn dirt can be removed with all commonly available aqueous solutions, naphtha or alcohol. For cleaning coated parts, use a linen or leather cloth that is moistened with one of these substances.

Caution:

Acetone, xylene or nitro-containing thinner can harm the microscope and thus must not be used.

Test cleaning solutions of unknown composition first on a less visible area of the unit. Be sure that coated or plastic surfaces do not become matted or etched.

Cleaning the stage

Rub the stage with paraffin oil or acid-free Vaseline to remove light spots on the stage.

10. Care of the microscope

Cleaning glass surfaces

Remove dust on glass surfaces with a fine, dry and fat-free hair brush, by blowing with a blow bag or vacuum suction.

Remove stubborn dirt on glass surfaces with a clean cloth dampened with distilled water. If the dirt still can not be removed, use pure alcohol, chloroform or benzene.

Cleaning objectives

Caution!

The objectives must not be unscrewed during cleaning. If damage appears on inner surfaces, the objectives must be sent to your Leica subsidiary for repair. We also advise against cleaning the inside surfaces of the eyepieces.

The front lenses of objectives are cleaned as described under "Cleaning Glass Surfaces". The upper lens is cleaned by being blown off with a pneumatic pump.

Removing immersion oil

Caution!

Follow safety instructions for immersion oil!

First, wipe off the immersion oil with a clean cotton cloth, and then re-wipe the surface several times with ethyl alcohol.

10.3 Handling Acids and Bases

For examinations using acids or other aggressive chemicals, particular care must be taken.

Caution:

Be absolutely certain to prevent the optics and mechanical parts from coming into contact with these chemicals.


11. Essential Wear and Spare Parts

Order No. Material No.	Name	Used for
<u>Replacement Lamp</u>		
11 500 974	Halogen lamp 12 V 100 W	107/2, 106, 107 lamp housing
11 500 138	High-pressure mercury burner 100 W	106 z lamp housing
11 500 321	High-pressure mercury burner 100 W (103 W/2)	106 z lamp housing
11 500 139	High-pressure xenon burner 75 W	106z lamp housing
<u>Screw cap for unused objective receptacles</u>		
020-422-570-000	Screw cap M 25	Objective turret
11 020-422-580-028	Screw cap M 32	Objective turret
<u>Cover for unused objective DIC disk opening</u>		
11 090-144-020-088	Cover for DIC	microscope stand
<u>Dust and light protection cover for analyzer slot</u>		
11 020-437-101-013	Analyzer slot cover	microscope stand
<u>Dust and light protection cover for camera port openings</u>		
11 020-387-556-009	Port slot cover	microscope stand
<u>Replacement eyecup (antiglare protection) for HC PLAN eyepiece</u>		
021-500-017-005	HC PLAN eyecup	10x/25 eyepiece
021-264-520-018	HC PLAN eyecup	10x/22 eyepiece
021-264-520-018	HC PLAN eyecup	10x/20 eyepiece
<u>Immersion oil conforming to DIN/ISO standards, fluorescence-free</u>		
11 513 787	10 ml	OIL and IMM objectives and oil condenser heads
11 513 522	100 ml	
11 513 788	500 ml	

12. Dimensions

12. Dimensions

Space requirements

Height compensation plate

A height compensation plate was developed to raise the viewing height by 20 mm or to raise the side camera ports for oversize cameras or spinning disks, or to use the microscope with an inactive bottom port on workbenches without openings. Two versions (12 mm and 25 mm) of the height compensation plate are available.

13. Abbreviations and Pictograms

	Contrasting method
	Magnification
	Illumination
	Ports/Eyepiece
	Focus
	Lower focus stop not set
	Lower focus stop set
	Focus position not set
	Focus position set
	Shutter open
	Shutter closed
	Transmitted-light filter
	Field diaphragm, rectangular
	Field diaphragm, round
	Aperture diaphragm
	Light distribution

13. Abbreviations and Pictograms

AP	Aperture diaphragm
BF	Bright field
CUBE	Reflector cube
DF	Dark field incident/transmitted light
DIC	Differential Interference Contrast
FD	Field diaphragm
FF	Focus Finder
FLUO	Fluorescence axis (incident light)
ICR	Interference contrast, incident light
ICT	Interference contrast, transmitted light
IL	Incident light
INT	Intensity
PH	Phase contrast
POL	Polarization, incident/transmitted light
RL	Incident light
TL	Transmitted light

14. Index

- Abbreviations** 95
- Active ports** 56
- Allen key** 20
- Ambient temperature** 9
- Analyzer** 35
- Analyzer IC/P** 33, 35
- Analyzer slot** 17
- Aperture diaphragm**
 - 14, 16, 57, 66, 86
- Aperture diaphragm adjustment** 61
- Arc** 64
- Assembly tools** 20
- Attachable mechanical stage** 17, 36
- Beam splitting** 73
- Bright field** 66, 88
- Camera** 48
- Centering keys** 20
- Centering the mercury lamps** 63
- Centering window** 17
- Cleaning glass surfaces** 92
- Cleaning objectives** 92
- Cleaning the stage** 91
- C-Mount 0.5x/0.63x** 48
- Coarse** 72
- Collector** 27
- Color coding (objectives)** 83
- Combined methods** 68
- Computer** 55
- Computer port** 50
- Condenser** 12, 14, 44, 46, 55
- Condenser base S1-S28** 44
- Condenser head S1** 44
- Condenser head S28** 44
- Condenser heads** 47
- Condenser height adjuster** 15
- Condenser prisms** 29, 44
- Condenser prisms** 44
- Condenser tool** 45
- Consumable parts** 93
- Contrast methods** 11, 56, 66
- Controls** 13, 51
- Correction for vision problems** 81
- Corrective mount** 83
- CUBE** 59
- Dark field** 67, 88
- Diaphragms** 56
- DIC module** 28
- DIC objective prism disk** 16
- DIC objective prisms** 28
- DIC specimen** 69
- Differential interference contrast** 69
- Digital camera** 48
- DIN VDE** 8
- Direct interface** 49
- Drawer** 17
- DRY** 82
- Dust cover** 91
- E-focus control buttons** 17
- Electromagnetic compatibility** 8
- Electronics box**
 - 9, 13, 15, 49, 55, 65
- EU directive** 8
- EXT1-EXT4 sockets** 40
- Eyebase** 73
- Eye piece identification** 74
- Eye piece tube** 16
- Eye pieces** 14, 16, 47, 74, 81
- Fast** 84
- Field diaphragm** 15, 86
- Field diaphragm**
 - 16, 57, 61, 66, 86
- Field diaphragm adjustment** 61
- Field diaphragm centering** 17
- Filter block** 32
- Filter cube** 32, 59
- Filter drawer** 59
- Filter slide** 16
- Filters** 15, 41
- Fine** 72
- Fixed function buttons**
 - 53, 57, 86
- Fixed stage** 36, 37
- Fluo drawer** 31
- Fluorescence** 70, 89
- Fluorescence filter cubes** 70
- Focus** 56
- Focus control buttons** 59, 71
- Focus Finder** 86

14. Index

- Focus position 71
Focus stop 59
Focus wheel 14, 16, 17
Focusing 71
Focusing telescope 39
Frequency 9
Front control panel 16, 59, 85
Function button 54
Function button assignment 54, 58
Function buttons 55, 57, 58, 59
Fuses 9
- G**
Graticule 78, 81
- H**
Halogen lamp 23, 42
Handling acids and bases 92
Height compensation plate 94
High-pressure mercury burner 100 W 25
High-pressure xenon burner 75 W 25
- I**
Illumination 87
Illumination telescope 21
IMM 82
Immersion objectives 82, 83
Immersion oil 93
Incident light lamp housing 17
Incident light polarizer L ICR/P 34
Incident light polarizers in slider 33
Incident Light Turret Disk 31
Incident light, fluorescence 61
Incident-light analyzer 35
Incident-light axis 11
Incident-light lamps 23
Incident-light polarizer R/P 33
- Initialization 65
Insert frame (coverslips) 39
Inserting discharge lamps 25
Inserting the lamp 43
Installation site 19
Intelligent automation 51
Intensity control 85
Intensity controls 14
Interfaces 13
- K**
Knobs 53, 60
- L**
Lamp bases 26
Lamp housing 106 z 25, 27, 62
Lamp housing 106 z L 24, 25
Lamp housing 107 42
Lamp housing 107 L 23, 62
Lamp housing 107/2 42, 62
Lamp housing for incident light 15
Lamp housing for transmitted light 15
Lamp mount 16, 21
Lamp replacement 41
LeicaScreen 16, 52, 55, 56, 65, 90
Light intensity 16, 57
Light rings 55
Light source adjustment 62
Light sources 85
Lower stop 59, 71
- M**
Magnification 56
Magnification changer 11, 59, 85
Manual method 68
Mechanical 3-plate stage 36
Medical instrument 8
Mirror housing 21
- Motorized 3 -plate or scanning stages 38
Motorized method 68
- O**
Object displacement 84
Objective 15, 56, 81
Objective changing 81
Objective identification 74
Objective turret 11, 15, 17, 30
Observation ports 12
Opener for drawer 17
Operating mode 82
Operating temperature 10
Overvoltage category 9
- P**
Parfocality compensation 30
PCI card (PC) 48
Phase contrast 70, 88
Phase rings 55
Phillips screwdriver 20
Pictograms 95
Polarization 68, 89
Pollution degree 9
Port 73
Port selection 73
Port switching 16
Power input: 9
Power supply 50
Precise 84
Protective gloves 21
Protective mask 21
- R**
Reflector BF 32
Reflector cubes 62
Reflector DF 32
Reflector ICR 32
Reflector Smith 32
Relative humidity: 9
Remote control module 17, 60

- Removing Immersion Oil 92
- Replacement eyecup 93
- Replacement lamps 93
- Replacement parts 93
- Reset function 52
- Rotating stage 39
- RS232 ports 49

- S**afety class 9
- Safety regulations 9
- Scanning stages 38
- Screw lengths 36
- Shutter 11, 59, 70
- Side-Port 16, 17
- SmartMove 17, 51, 53
- Software 53
- Software tools 13
- Specifications 9
- Specimen stage 15, 36
- Stage positions 84
- Stages 12, 17, 84
- Stand 14, 87
- Stand package 18
- Step increments 72, 84
- Stops 56, 71
- Supply unit ebq 100 9, 27, 50
- Supply voltage 9
- Switching on the microscope 55
- System package 18


- T**L/IL switching 16
- Toggling transmitted light/incident light 57
- Top port 16
- Transmitted-light axis 11
- Transmitted-light illumination carrier 40
- Transmitted-light illuminator 15
- Transmitted-light lamp housing 41
- Transmitted-light specimen 66, 70
- Transportation 19
- Tube 11, 14
- Tube setting 73
- Turret disk 32

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Copyright © Leica Microsystems Wetzlar GmbH, Ernst-Leitz-Straße, 35578 Wetzlar, Germany 2002. Tel. (064 41) 29-0, Fax (064 41) 29-25 99. LEICA and the Leica logos are registered trademarks of Leica IR GmbH.